



# **Aqua-Duct Documentation**

*Release 0.5.15*

**Tomasz Magdziarz      Karolina Mitusińska**  
**Agata Raczyńska      Artur Góra**

**Jul 27, 2019**



# CONTENTS

<b>1</b>	<b>Aqua-Duct installation guide</b>	<b>1</b>
1.1	Overview	1
1.2	Troubleshooting	1
1.3	Requirements	1
1.3.1	Software-wise requirements	1
1.3.2	Hardware-wise requirements	2
1.4	Installation	2
1.4.1	Generic Python installation	2
1.4.2	GNU/Linux	2
1.4.3	macOS	3
1.4.4	Windows	4
1.4.5	OpenBSD	5
<b>2</b>	<b>Valve manual</b>	<b>7</b>
2.1	Valve invocation	7
2.1.1	Usage	7
2.1.2	Configuration file template	8
2.1.3	Valve calculation run	8
2.2	How does Valve work	9
2.2.1	Traceable residues	9
2.2.2	Raw paths	11
2.2.3	Separate paths	11
2.2.4	Clusterization of inlets	12
2.2.5	Passing paths	14
2.2.6	Analysis	14
2.2.7	Visualization	19
<b>3</b>	<b>Configuration file options</b>	<b>23</b>
3.1	Section <b>global</b>	23
3.2	Common settings of stage sections	24
3.3	Stage <b>traceable_residues</b>	24
3.4	Stage <b>raw_paths</b>	25
3.5	Stage <b>separate_paths</b>	26
3.6	Stage <b>inlets_clusterization</b>	27
3.7	Stage <b>analysis</b>	27
3.8	Stage <b>visualize</b>	27
3.9	Clusterization sections	29
3.9.1	barber	30
3.9.2	dbscan	30
3.9.3	affprop	31
3.9.4	meanshift	31
3.9.5	birch	32
3.9.6	kmeans	32
3.10	Smooth section	32

<b>4</b>	<b>Valve tutorial</b>	<b>33</b>
4.1	Valve invocation . . . . .	33
4.2	Test data . . . . .	33
4.3	Inspect your system . . . . .	33
4.3.1	Create <i>Object definition</i> . . . . .	33
4.3.2	Create <i>Scope definition</i> . . . . .	34
4.4	Prepare config file . . . . .	34
4.5	Run <i>Valve</i> . . . . .	34
4.5.1	Visual inspection . . . . .	34
4.5.2	Clusterization . . . . .	34
4.5.3	Analysis tables . . . . .	34
4.6	Feedback . . . . .	34
<b>5</b>	<b>aquaduct</b>	<b>37</b>
5.1	aquaduct package . . . . .	37
5.1.1	Subpackages . . . . .	37
5.1.2	Module contents . . . . .	74
<b>6</b>	<b>Aqua-Duct changelog</b>	<b>75</b>
	<b>Python Module Index</b>	<b>77</b>
	<b>Index</b>	<b>79</b>

## AQUA-DUCT INSTALLATION GUIDE

### 1.1 Overview

Aqua-Duct software is software written in Python (CPython) and comprises of two elements:

1. aquaduct - a Python package,
2. valve - a script that uses *aquaduct* to perform calculations.

#### Download

You can download Aqua-Duct packages directly from [Aqua-Duct homepage](#). This page includes older versions of Aqua-Duct as well as development version.

If you follow this installation guide you will install current release.

### 1.2 Troubleshooting

If you encounter any problems with installation do not hesitate to contact us at [info@aquaduct.pl](mailto:info@aquaduct.pl). We are **RE-ALLY** willing to help!

Please, provide us with as much info as you can. In particular try to include following information:

- Operating system's name and version, and CPU architecture (if relevant).
- Python version.
- Command(s) you have used for installation.
- Any error/warning/info message(s) that emerged during or after installation.

### 1.3 Requirements

#### 1.3.1 Software-wise requirements

- numpy  $\geq$  1.10.0
- scipy  $\geq$  0.17.1
- scikit-learn  $\geq$  0.16.0
- MDAnalysis[amber]  $\geq$  0.16.0,  $<$  0.17.0
- joblib  $\geq$  0.10

## 1.3.2 Hardware-wise requirements

Aqua-Duct should work on every machine on which you can install the above mentioned software. On computers older than 10 years it may work very slow though. We recommend 64bit SMP architecture, with at least 4GB RAM (32 GB RAM is recommended).

## 1.4 Installation

### 1.4.1 Generic Python installation

The easiest way to install Aqua-Duct is to install Python 2.7 and use following command:

```
pip install aquaduct
```

If *pip* is not available try to install it by typing:

```
easy_install pip
```

Depending on the settings of your system you can prepend the above command with *sudo* or *doas* or do *user* installation:

```
# sudo
sudo pip install aquaduct

# doas
doas pip install aquaduct

# 'user' installation
pip install aquaduct --user
```

It is also good idea to try to install Aqua-Duct using virtualenv:

```
virtualenv aquaduct_installation
cd aquaduct_installation
. bin/activate
pip install aquaduct
```

#### 1.4.1.1 Installation of PyMOL

Under most modern GNU/Linux distributions PyMOL is available as a package in repositories. For example if you are under Ubuntu/Debian you can install it by following command:

```
sudo apt-get install pymol
```

Under Windows there are several ways to install PyMOL, for more details see [PyMOL web site](#).

Instructions for macOS and OpenBSD are in appropriate sections below.

### 1.4.2 GNU/Linux

Installation was tested on limited number of GNU/Linux systems. On the most of modern installations you can simply follow generic instructions, for example under Ubuntu 16.04 you can type:

```
sudo pip install aquaduct
```

### 1.4.2.1 NetCDF4 & MDAAnalysis installation Ubuntu 14.04

Other systems may require additional work, in particular installation of NetCDF4 is sometimes cumbersome. Following is an example how to install all required packages under Ubuntu 14.04:

```
# install required python packages
sudo apt-get install python-dev python-pip python-numpy python-scipy python-
↳matplotlib python-scikits-learn

# install necessary libraries and git - all required to compile netCDF4
sudo apt-get -y install libnetcdf-dev libhdf5-dev git

# clone netcdf4 python repository
git clone https://github.com/Unidata/netcdf4-python.git
# cd to cloned repository
cd netcdf4-python
# modify setup.cfg to add paths of hdf5 and netcdf4 libraries
sed -i '/\[directories\]/a \
HDF5_dir = /usr/lib \
HDF5_libdir = /usr/lib \
HDF5_incl_dir = /usr/include \
netCDF4_dir = /usr/lib \
netCDF4_libdir = /usr/lib \
netCDF4_incl_dir = /usr/include' setup.cfg
# run setup.py
sudo python setup.py install

# install MDAAnalysis
sudo pip install "MDAAnalysis[amber]==0.16.2"
```

If everything went fine you can follow generic instructions.

### 1.4.2.2 SciPy update and Ubuntu/Debian

Debian (and Ubuntu) uses strange approach to Python installation. To install newer version of SciPy (if required) try following procedure:

```
# install libraries required for SciPy compilation
apt-get build-dep python-scipy

# install SciPy
easy_install-2.7 --upgrade scipy
```

**Warning:** The above procedure will remove current SciPy from *easy-install.pth* file.

## 1.4.3 macOS

Aqua-Duct installation was tested on macOS Sierra and is quite straightforward. It can be installed either with existing system Python or with custom Python installation. In both cases one have to install Xcode for the App Store.

### 1.4.3.1 System native Python

```
sudo easy_install pip
sudo pip install aquaduct
```

The drawback of using system Python installation is a lack of PyMOL. It should be, however, relatively easy to compile PyMOL on your own. Try to follow compilation instruction under BSD systems.

### 1.4.3.2 Custom Python

This is recommended way of Aqua-Duct installation. If you do not have custom Python installation you can get it by using one of package managers available for macOS, for example [homebrew](#). With this package manager you can do following:

```
brew install python
sudo easy_install pip
sudo pip install aquaduct
```

Next, you can install PyMOL:

```
brew install pymol
brew cask install xquartz
```

Once XQuartz is installed you should reboot. The above procedure installs PyMOL, however, PyMOL Python modules are not visible. To fix it you can issue following commands:

```
cd /usr/local/lib/python2.7/site-packages
sudo ln -s /usr/local/Cellar/pymol/*/libexec/lib/python2.7/site-packages/* ./
```

The above instruction assumes that you are using *brew* and you have only one PyMOL installation.

## 1.4.4 Windows

Installation under Windows is also possible. The limiting factor is MDAnalysis which is not officially available under Windows yet. You can, however, install Cygwin and perform Aqua-Duct installation in Cygwin.

First, start with [Cygwin installation](#). During the setup select following packages:

- python (2.7)
- python-devel (2.7)
- python-cython
- libnetcdf-devel
- libhdf5-devel
- liblapack-devel
- libopenblas
- python-numpy
- python-six

Another key component that have to be installed is C, C++ and Fortran compilers. You can simply install **gcc-g++** and **gcc-fortran** packages as a first choice, select following packages:

- gcc-g++
- gcc-fortran

Once Cygwin is installed with all required libraries you can perform following steps:

```
# install pip
easy_install-2.7 pip
```

First, try to install SciPy:



```
# install SciPy
pip install scipy
```

If you encounter any problems related to missing **xlocale.h** header file try the following workaround:

```
# prepare fake xlocale.h
ln -s /usr/include/locale.h xlocale.h
export CFLAGS="I"$( pwd )

# install SciPy
pip install scipy
```

---

**Note:** The above procedure for SciPy installation might not be optimal. For more information please go to [SciPy web page](#).

---

Now, install **scikit-learn** and then Aqua-Duct:

```
# install scikit-learn
pip install scikit-learn

# finally, install aquaduct
pip install aquaduct
```

## 1.4.5 OpenBSD

Aqua-Duct can be also installed under OpenBSD (5.9 and 6.0 amd64). NetCDF-c version 4 has to be installed as OpenBSD ships only netCDF in version 3. First, install hdf5 library and GNU make:

```
# install hdf5 and GNU make
pkg_add hdf5 gmake
```

Next, download netCDF sources. Version 4.2.1.1 works out of the box but is a bit outdated. Visit [NetCDF web page](#) and select version of your choice. Older versions are available in the [FTP archive](#). Once netCDF is downloaded and extracted go to the source directory and try following procedure:

```
# set LD and CPP flags
export LDFLAGS=-L/usr/local/lib
export CPPFLAGS=-I/usr/local/include

# configure project
./configure --enable-shared --enable-dap --disable-doxygen --enable-netcdf-4 --
↪prefix=/path/to/netCDF4/lib

# make and install
gmake
gmake install
```

You may now install py-scipy package:

```
pkg_add py-scipy
```

Install pip if it is missing:

```
pkg_add py-pip
```

Install netCDF4 Python:

```
# define netcdf-4 installation directory
export NETCDF4_DIR=/path/to/netCDF4/lib
pip2.7 install netCDF4
```

At this point you can follow generic Python instructions, type:

```
pip2.7 install aquaduct
```

### 1.4.5.1 PyMOL at OpenBSD

According to our knowledge it is possible to install PyMOL 1.4.1 and it is sufficient to work with Aqua-Duct. Go to [SourceForge PyMOL download page](#) and download, save, and extract sources.

PyMOL requires Python Mega Widgets. Download, for example Pmw 1.3.3b from [SourceForge Pmw download page](#). Extract it and install by:

```
python2.7 setup.py install
```

TKinter (2.7) and several other packages are also required:

```
pkg_add python-tkinter freeglut glew png
```

Next, go to the extracted PyMOL sources open setup.py and modify inc\_dirs variable at line 129 by adding following paths:

```
"/usr/X11R6/include/freetype2",
"/usr/X11R6/include",
"/usr/local/include",
```

Now, you can build and install PyMOL by typing following commands:

```
python2.7 setup.py build
python2.7 setup.py install
python2.7 setup2.py install
cp pymol /usr/local/bin
```

PyMOL can be run by typing *pymol* or can be used as Python module.

### 1.4.5.2 Other BSDs

Installation on other BSDs might be easier. For example, Python netCDF4 is available in ports of FreeBSD and DragonFlyBSD. Try to install it and SciPy, then proceed to generic Python installation instructions.

If you are using NetBSD or other BSD try to follow OpenBSD instructions.

## VALVE MANUAL

Valve application is a driver that uses *aquaduct* module to perform analysis of trajectories of selected residues in Molecular Dynamics simulation.

### 2.1 Valve invocation

Once *aquaduct* module is installed (see *Aqua-Duct installation guide*) properly on the machine, Valve is available as `valve.py` command line tool.

#### 2.1.1 Usage

Basic help of Valve usage can be displayed by following command:

```
valve.py --help
```

It should display following information:

```
usage: valve.py [-h] [--debug] [--debug-file DEBUG_FILE]
               [--dump-template-config] [-t THREADS] [-c CONFIG_FILE] [--sps]
               [--max-frame MAX_FRAME] [--min-frame MIN_FRAME]
               [--step-frame STEP_FRAME] [--sandwich] [--cache-dir CACHEDIR]
               [--cache-mem] [--version] [--license]
```

Valve, Aquaduct driver

optional arguments:

<code>-h, --help</code>	show this help message <b>and</b> exit
<code>--debug</code>	Prints debug info. (default: <b>False</b> )
<code>--debug-file DEBUG_FILE</code>	Debug log file. (default: <b>None</b> )
<code>--dump-template-config</code>	Dumps template config file. Suppress <b>all</b> other output <b>or</b> actions. (default: <b>False</b> )
<code>-t THREADS</code>	Limit Aqua-Duct calculations to given number of threads. (default: <b>None</b> )
<code>-c CONFIG_FILE</code>	Config file filename. (default: <b>None</b> )
<code>--sps</code>	Use single precision to store data. (default: <b>False</b> )
<code>--max-frame MAX_FRAME</code>	Maximal number of frame. (default: <b>None</b> )
<code>--min-frame MIN_FRAME</code>	Minimal number of frame. (default: <b>None</b> )
<code>--step-frame STEP_FRAME</code>	Frames step. (default: <b>None</b> )
<code>--sandwich</code>	Sandwich mode <b>for</b> multiple trajectories. (default: <b>False</b> )
<code>--cache-dir CACHEDIR</code>	Directory <b>for</b> coordinates caching. (default: <b>None</b> )

(continues on next page)

(continued from previous page)

<code>--cache-mem</code>	Switch on memory caching. (default: <b>False</b> )
<code>--version</code>	Prints versions <b>and</b> exits. (default: <b>False</b> )
<code>--license</code>	Prints short license info <b>and</b> exits. (default: <b>False</b> )

## 2.1.2 Configuration file template

Configuration file used by *Valve* is of moderate length and complexity. It can be easily prepared with a template file that can be printed by *Valve*. Use following command to print configuration file template on the screen:

```
valve.py --dump-template-config
```

Configuration file template can also be easily saved in to a file with:

```
valve.py --dump-template-config > config.txt
```

Where config.txt is a configuration file template.

For detailed description of configuration file and available options see [Configuration file options](#).

## 2.1.3 Valve calculation run

Once configuration file is ready *Valve* calculations can be run with a following simple command:

```
valve.py -c config.txt
```

Some of *Valve* calculations can be run in parallel. By default all available CPU cores are used. This is not always desired - limitation of used CPU cores can be done with `-t` option which limits number of concurrent threads used by *Valve*. If it equals 1 no parallelism is used.

---

**Note:** Specifying number of threads greater then available CPU cores is generally not optimal.

However, in order to maximize usage of available CPU power it is recommended to set it as number of cores + 1. The reason is that *Valve* uses one thread for the main process and the excess over one for processes for parallel calculations. When parallel calculations are executed the main thread waits for results.

---

---

**Note:** Options `--min-frame`, `--max-frame`, and `--step-frame` can be used to limit calculations to specific part of trajectory. For example, to run calculations for 1000 frames starting from frame 5000 use following options: `--min-frame 4999 --max-frame 5999`; to run calculations for every 5th frame use: `--step-frame 5`.

---

### 2.1.3.1 Single precision storage

Most of the calculation is *Valve* is performed by NumPy. By default, NumPy uses double precision floats. *Valve* does not change this behavior but has special option `--sps` which forces to store all data (both internal data stored in RAM and on the disk) in single precision. This spare a lot of RAM and is recommended what you perform calculation for long trajectories and you have limited amount of RAM.

### 2.1.3.2 Cache

Storage of coordinates for all paths for very long MD trajectories requires huge amount of RAM. User can decide whether *aqueduct* should store coordinates in memory or in separated directory. Option `--cache-mem`

instruct *Valve* to store coordinates in RAM; `--cache-dir` stores coordinates in selected directory. If neither of both options is selected, coordinates are calculated on demand.

---

**Note:** If no cache is used (memory or dir) *Master paths* cannot be calculated.

---

### 2.1.3.3 Sandwich

Trajectory data can be provided as several files. By default these files are processed in sequential manner making one long trajectory. If option `--sandwich` is used trajectory files are read as layers. For each layer, search of traceable residues is done separately (stage I and II) but processing and analysis (stage III, IV, V, and VI) are done for all paths simultaneously. Usage of `--sandwich` option is further referenced as *sandwich* mode.

### 2.1.3.4 Debugging

*Valve* can output some debug information. Use `--debug` to see all debug information on the screen or use `--debug-file` with some file name to dump all debug messages to the given file. Beside debug messages standard messages will be saved in the file as well.

## 2.2 How does *Valve* work

Application starts with parsing input options. If `--help` or `--dump-template-config` options are used appropriate messages are printed on the screen and *Valve* quits with signal 0.

---

**Note:** In current version *Valve* does not check the validity of the config file.

---

If config file is provided (option `-c`) *Valve* parse it quickly and regular calculations starts according to its content. Calculations performed by *Valve* are done in six stages described in the next sections.

### 2.2.1 Traceable residues

In the first stage of calculation *Valve* finds all residues that should be traced and appends them to the list of *traceable residues*. It is done in a loop over all frames. In each frame residues of interest are searched and appended to the list but only if they are not already present on the list. In *sandwich\_option* mode this is repeated for each layer.

The search of *traceable residues* is done according to user provided specifications. Two requirements have to be met to append residue to the list:

1. The residue has to be found according to the *object* definition.
2. The residue has to be within the *scope* of interest.

The *object* definition encompasses usually the active site of the protein (or other region of interest of macro-molecule in question). The *scope* of interest defines, on the other hand, the boundaries in which residues are traced and is usually defined as protein.

Since *aqueduct* in its current version uses *MDAnalysis* Python module for reading, parsing and searching of MD trajectory data, definitions of *object* and *scope* have to be given as its *Selection Commands*.

#### 2.2.1.1 Object definition

*Object* definition has to comprise of two elements:

1. It has to define residues to trace.

2. It has to define spatial boundaries of the *object* site.

For example, proper *object* definition could be following:

```
(resname WAT) and (sphzone 6.0 (resnum 99 or resnum 147))
```

It defines WAT as residues that should be traced and defines spatial constraints of the *object* site as spherical zone within 6 Angstroms of the center of masses of residues with number 99 and 147.

### 2.2.1.2 Scope definition

*Scope* can be defined in two ways: as *object* but with broader boundaries or as the convex hull of selected molecular object.

In the first case definition is very similar to *object* and it has to follow the same limitations. For example, proper *scope* definition could be following:

```
resname WAT and around 2.0 protein
```


It consequently has to define WAT as residues of interest and defines spatial constraints: all WAT residues that are within 2 Angstroms of the protein.

If the *scope* is defined as the convex hull of selected molecular object (which is **recommended**), the definition itself have to comprise of this molecular object only, for example *protein*. In that case the scope is interpreted as the interior of the convex hull of atoms from the definition. Therefore, *traceable residues* would be in the scope only if they are within the convex hull of atoms of *protein*.

### Convex hulls of macromolecule atoms

AQ uses quickhull algorithm for convex hulls calculations (via SciPy class `scipy.spatial.ConvexHull`, see also <http://www.qhull.org/> and original publication [The quickhull algorithm for convex hulls](#)).

Convex hull concept is used to check if traced molecules are inside of the macromolecule. Convex hull can be considered as rough approximation of molecular surface. Following picture shows schematic comparison of convex hull and solvent excluded surface (SES):



Convex hull (red shape) of atoms (blue dots with VdW spheres) and SES (blue line): a) convex hull and SES cover roughly the same area, convex hull approximates SES; b) movement of one atom dramatically changes SES, however, interior of the molecule as approximated by convex hull remains stable.

No doubts, convex hull is a very rough approximation of SES. It has, however, one very important property when it is used to approximate interior of molecules: its interior does not considerably depend on the molecular conformation of a molecule (or molecular entity) in question.

## 2.2.2 Raw paths

The second stage of calculations uses the list of all traceable residues from the first stage and for each residue in each frame two checks are performed:

1. Is the residue in the *scope* (this is always calculated according to the scope definition).
2. Is the residue in the *object*. This information is partially calculated in the first stage and can be reused in the second. However, it is also possible to recalculate this data according to the new *object* definition.

For each of the *traceable residues* a special *Path* object is created which stores frames in which a residue is in *scope* or in *object*.

---

**Note:** Residue is in *object* only if it is also in *scope*.

---

## 2.2.3 Separate paths

The third stage uses collection of *Path* objects to create *Separate Path* objects. Each *Path* comprise data for one residue. It may happen that the residue enters and leaves the *scope* and the *object* many times over the entire MD. Each such event is considered by *Valve* as a separate path.

There are two types of *Separate Paths*:

- *Object Paths*
- *Passing Paths*

*Object Paths* are traces of molecules that visited *Object* area. *Passing Paths* are traces of molecules that entered *Scope* but did not entered *Object* area.

*Passing paths* comprises of one part only. Each *object path* comprises of three parts:

1. *Incoming* - Defined as a path that leads from the point in which residue enters the *scope* and enters the *object* for the first time.
2. *Object* - Defined as a path that leads from the point in which residue enters the *object* for the first time and leaves it for the last time.
3. *Outgoing* - Defined as a path that leads from the point in which residue leaves the *object* for the last time and leaves the *scope*.

It is also possible that incoming and/or outgoing part of the separate path is empty.

---

**Note:** Generation of *Passing paths* is optional and can be switched off.

---

**Warning:** Generation of *Passing paths* without redefinition of *Object* area in stage I and II may lead to false results.

### 2.2.3.1 Auto Barber

After the initial search of *Separate Path* objects it is possible to run Auto Barber procedure which trims paths down to the approximated surface of the macromolecule or other molecular entity defined by the user. This trimming is done by creating collection of spheres that have centers at the ends of paths and radii equal to the distance for the center to the nearest atom of user defined molecular entity. Next, parts of raw paths that are inside these spheres are removed and separate paths are recreated.

Auto Barber procedure has several options, for example:

- **auto\_barber** allows to define molecular entity which is used to calculate radii of spheres used for trimming raw paths.
- **auto\_barber\_mincut** allows to define minimal radius of spheres. Spheres of radius smaller then this value are not used in trimming.
- **auto\_barber\_maxcut** allows to define maximal radius of spheres. Spheres of radius greater then this value are not used in trimming.
- **auto\_barber\_tovdw** if set to *True* radii of spheres are corrected (decreased) by Van der Waals radius of the closest atom.

See also *options of separate\_paths* stage.

### 2.2.3.2 Smoothing

Separate paths can be optionally smoothed. Current *aqueduct* version allows perform *soft* smoothing only, i.e. smoothing is used only for visualization purposes. Raw paths cannot be replaced by the smoothed.

### Available methods

Aqua-Duct implements several smoothing methods:

1. Savitzky-Golay filter - *SavgolSmooth* - see also original publication [Smoothing and Differentiation of Data by Simplified Least Squares Procedures](#) (doi:10.1021/ac60214a047).
2. Window smoothing - *WindowSmooth*
3. Distance Window smoothing - *DistanceWindowSmooth*
4. Active Window smoothing - *ActiveWindowSmooth*
5. Max Step smoothing - *MaxStepSmooth*
6. Window over Max Step smoothing - *WindowOverMaxStepSmooth*
7. Distance Window over Max Step smoothing - *DistanceWindowOverMaxStepSmooth*
8. Active Window over Max Step smoothing - *ActiveWindowOverMaxStepSmooth*

For detailed information on available configuration options see configuration file *smooth section* description.

### 2.2.4 Clusterization of inlets

Each of the separate paths has beginning and end. If they are at the boundaries of the *scope* they are considered as *Inlets*, i.e. points that mark where the *traceable residues* enters or leaves the *scope*. Clusters of inlets, on the other hand, mark endings of tunnels or ways in the system which was simulated in the MD.

Clusterization of inlets is performed in following steps:

1. *Initial clusterization*: All inlets are submitted to selected clusterization method and depending on the method and settings, some of the inlets might not be arranged to any cluster and are considered as outliers.
2. [Optional] *Outliers detection*: Arrangement of inlets to clusters is sometimes far from optimal. In this step, *inlets* that do not fit to cluster are detected and annotated as outliers. This step can be executed in two modes:
  1. *Automatic mode*: Inlet is considered to be an outlier if its distance from the centroid is greater than mean distance + 4 \* standard deviation of all distances within the cluster.
  2. *Defined threshold*: Inlet is considered to be an outlier if its minimal distance from any other point in the cluster is greater than the threshold.



3. [Optional] *Reclusterization of outliers*: It may happen that the outliers form actually clusters but it was not recognized in initial clusterization. In this step clusterization is executed for outliers only and found clusters are appended to the clusters identified in the first step. Rest of the inlets are marked as outliers.

#### 2.2.4.1 Potentially recursive clusterization

Both *Initial clusterization* and *Reclusterization* can be run in a recursive manner. If in the appropriate sections defining clusterization methods option *recursive\_clusterization* is used appropriate method is run for each cluster separately. Clusters of specific size can be excluded from recursive clusterization (option *recursive\_threshold*). It is also possible to limit maximal number of recursive levels - option *max\_level*.

For additional information see *clusterization sections* options.

#### 2.2.4.2 Available methods

Aqua-Duct implements several clustering methods. The recommended method is **barber** method which bases on *Auto Barber* procedure. Rest of the methods are implemented with `sklearn.cluster` module:

1. `aquaduct.geom.cluster.BarberCluster` - default for *Initial clusterization*. It gives excellent results. For more information see *barber clusterization method* description.
2. `MeanShift` - see also original publication *Mean shift: a robust approach toward feature space analysis* (doi:10.1109/34.1000236).
3. `DBSCAN` - default for *Reclusterization of outliers*, see also original publication *A Density-Based Algorithm for Discovering Clusters in Large Spatial Databases with Noise*
4. `AffinityPropagation` - see also original publication *Clustering by Passing Messages Between Data Points* (doi:10.1126/science.1136800)
5. `KMeans` - see also *k-means++: The advantages of careful seeding*, Arthur, David, and Sergei Vassilvitskii in *Proceedings of the eighteenth annual ACM-SIAM symposium on Discrete algorithms, Society for Industrial and Applied Mathematics* (2007), pages 1027-1035.
6. `Birch` - see also Tian Zhang, Raghu Ramakrishnan, Maron Livny *BIRCH: An efficient data clustering method for large databases* and Roberto Perdisci *JBirch - Java implementation of BIRCH clustering algorithm*.

For additional information see *clusterization sections* options.

#### 2.2.4.3 Master paths

At the end of clusterization stage it is possible to run procedure for *master path* generation. First, separate paths are grouped according to clusters. Paths that begin and end in particular clusters are grouped together. Next, for each group a *master path* (i.e., average path) is generated in following steps:

1. First, length of *master path* is determined. Lengths of each parts (incoming, object, outgoing) for each separate paths are normalized with bias towards longest paths. These normalized lengths are then used for as weights in averaging not normalized lengths. Values for all parts are summed and resulting value is the desired length of *master path*.
2. All separate paths are divided into chunks. Number of chunks is equal to the desired length of *master path* calculated in the previous step. Lengths of separate paths can be quite diverse, therefore, for different paths chunks are of different lengths.
3. For each chunk averaging procedure is run:
  1. Coordinates for all separate paths for given chunk are collected.
  2. Normalized lengths with bias toward longest paths for all separate paths for given chunk are collected.
  3. New coordinates are calculated as weighted average of collected coordinates. As weights collected normalized lengths are used.

4. In addition width of chunk is calculated as a mean value of collected coordinates mutual distances.
5. Type of chunk is calculated as probability (frequency) of being in the *scope*.
4. Results for all chunks are collected, types probability are changed to types. All data is then used to create Master Path. If this fails no path is created.

More technical details on master path generation can be found in `aquaduct.geom.master.CTypeSpathsCollection.get_master_path()` method documentation.

## 2.2.5 Passing paths

If *Passing paths* are allowed (see `allow_passing_paths` option in *separate\_paths configuration*) they will be generated using list of *traceable residues* from the first stage of calculations. In usual settings, where *Object* and *Scope* definitions are the same in both I and II stage, this will result in relatively low number of passing paths. In particular this will not show the real number of traced molecules that enter *Scope* during the simulation.

To get correct picture following options and settings have to be considered:

- **Stage `traceable_residues`**
  - `object` should be broad enough to encompass all molecules that should be traced. For example, if water is traced, `object` definition could be following: `resname WAT`.
- **Stage `raw_paths`**
  - In order to retain default Aqua-Duct behavior of tracing molecules that flow through *Object* area, it have to be redefined to encompass the active site only - see *Object definition* discussion.
  - `clear_in_object_info` should be set to `True`. Otherwise, traceable molecules will be limited according to current `object` definition but *Object* boundaries from **`traceable_residues`** stage will be used.
- **Stage `separate_paths`**
  - `allow_passing_paths` should be set to `True`. This allows generation of passing paths.

Additionally, in stage **`inlets_clusterization`** following options could also be adjusted:

- `exclude_passing_in_clusterization` could be set to `True`. This will exclude passing paths inlets from clusterization.
- If passing paths are not clustered they will be added as outliers. Option `add_passing_to_clusters` allows to add some of passing paths inlets to already existing clusters. This is done by Auto Barber method and therefore this option should define molecular entity used in Auto Barber procedure, for example `protein`.

## 2.2.6 Analysis

Fifth stage of *Valve* calculations analyses results calculated in stages 1 to 4. Results of the analysis are displayed on the screen or can be saved to text file and comprise of several parts.

### 2.2.6.1 General summary

Results starts with general summary.

- Title and data stamp.
- [Optional] Dump of configuration options.
- Frames window.

- **Names of traced molecules.**

---

**Note:** If more than one name is on the list all consecutive sections of *Analysis* results are provided for each name separately and, as well as, for all names.

---

- Number of traceable residues.
- Number of separate paths.
- Number of inlets.
- **Number of clusters.**
  - Outliers flag, *yes* if they are present.
- Clustering history - a tree summarizing calculated clusters.

### 2.2.6.2 Clusters statistics

- **Clusters summary - inlets.**
  - **Summary of inlets clusters. Table with 4 columns:**
    1. **Cluster:** ID of the cluster. Outliers have 0.
    2. **Size:** Size of the cluster.
    3. **INCOMING:** Number of inlets corresponding to separate paths that enter the scope.
    4. **OUTGOING:** Number of inlets corresponding to separate paths that leave the scope.
- **Cluster statistics.**
  - **Probabilities of transfers. Table with 7 columns:**
    1. **Cluster:** ID of the cluster. Outliers have 0.
    2. **IN-OUT:** Number of separate paths that both enter and leave the scope by this cluster.
    3. **diff: Number of separate paths that:**
      - \* Enter the scope by this cluster but leave the scope by another cluster, or
      - \* Enter the scope by another cluster but leave the scope by this cluster.
    4. **N: Number of separate paths that:**
      - \* Enter the scope by this cluster and stays in the object, or
      - \* Leaves the scope by this cluster after staying in the object.
    5. **IN-OUT\_prob:** Probability of **IN-OUT**.
    6. **diff\_prob:** Probability of **diff**.
    7. **N\_prob:** Probability of **N**.
  - **Mean lengths of transfers. Table with 8 columns:**
    1. **Cluster:** ID of the cluster. Outliers have 0.
    2. **X->Obj:** Mean length of separate paths leading from this cluster to the object.
    3. **Obj->X:** Mean length of separate paths leading from the object to this cluster.
    4. **p-value:** p-value of *ttest* of comparing **X->Obj** and **Obj->X**.
    5. **X->ObjMin:** Minimal value of length of separate paths leading from this cluster to the object.
    6. **X->ObjMinID:** ID of separate path for which **X->ObjMin** was calculated.

7. **Obj->XMin**: Minimal value of length of separate paths leading from the object to this cluster.
  8. **Obj->XMinID**: ID of separate path for which **Obj->XMin** was calculated.
- **Mean frames numbers of transfers. Table with 8 columns:**
1. **Cluster**: ID of the cluster. Outliers have 0.
  2. **X->Obj**: Mean number of frames of separate paths leading from this cluster to the object.
  3. **Obj->X**: Mean number of frames of separate paths leading from the object to this cluster.
  4. **p-value**: p-value of *ttest* of comparing **X->Obj** and **Obj->X**.
  5. **X->ObjMin**: Minimal value of number of frames of separate paths leading from this cluster to the object.
  6. **X->ObjMinID**: ID of separate path for which **X->ObjMin** was calculated.
  7. **Obj->XMin**: Minimal value of number of frames of separate paths leading from the object to this cluster.
  8. **Obj->XMinID**: ID of separate path for which **Obj->XMin** was calculated.

---

**Note:** Distributions of **X->Obj** and **Obj->X** might be not normal, *ttest* may result unrealistic values. This test will be changed in the future releases.

---

### 2.2.6.3 Clusters types statistics

- **Separate paths clusters types summary. Tables with 11 columns.**

– **Mean length of paths:**

1. **CType**: Separate path Cluster Type.
2. **Size**: Number of separate paths belonging to Cluster type.
3. **Size%**: Percentage of **Size** relative to the total number of separate paths.
4. **Tot**: Average total length of paths.
5. **TotStd**: Standard deviation of **Tot**.
6. **Inp**: Average length of incoming part of paths. If no incoming parts are available, NaN is printed (not a number).
7. **InpStd**: Standard deviation of **Inp**.
8. **Obj**: Average length of object part of paths. If no incoming parts are available, NaN is printed.
9. **ObjStd**: Standard deviation of **Inp**.
10. **Out**: Average length of outgoing part of paths. If no incoming parts are available, NaN is printed.
11. **OutStd**: Standard deviation of **Inp**.

– **Mean number of frames:**

1. **CType**: Separate path Cluster Type.
2. **Size**: Number of separate paths belonging to Cluster type.
3. **Size%**: Percentage of **Size** relative to the total number of separate paths.
4. **Tot**: Average total number of frames of paths.
5. **TotStd**: Standard deviation of **Tot**.

6. **Inp**: Average total number of incoming part of paths. If no incoming parts are available, NaN is printed (not a number).
7. **InpStd**: Standard deviation of **Inp**.
8. **Obj**: Average total number of object part of paths. If no incoming parts are available, NaN is printed.
9. **ObjStd**: Standard deviation of **Inp**.
10. **Out**: Average total number of outgoing part of paths. If no incoming parts are available, NaN is printed.
11. **OutStd**: Standard deviation of **Inp**.

### Cluster Type of separate path

Clusters types (or CType) is a mnemonic for separate paths that leads from one cluster to another, including paths that start/end in the same cluster or start/end in the *Object* area.

Each separate path has two ends: beginning and end. Both of them either belong to one of the clusters of inlets, or are among outliers, or are inside the scope. If an end belongs to one of the clusters (including outliers) it has ID of the cluster. If it is inside the scope it has special ID of N. Cluster type is an ID composed of IDs of both ends of separate path separated by colon charter.

#### 2.2.6.4 All separate paths data

- **List of separate paths and their properties. Table with 20 columns.**

1. **ID**: - Separate path ID.
2. **RES**: - Residue name.
3. **BeginF**: Number of frame in which the path begins.
4. **InpF**: Number of frame in which path begins Incoming part.
5. **ObjF**: Number of frame in which path begins Object part.
6. **OutF**: Number of frame in which path begins Outgoing part.
7. **EndF**: Number of frame in which the path ends.
8. **TotL**: Total length of path.
9. **InpL**: Length of Incoming part. If no incoming part NaN is given.
10. **ObjL**: Length of Object part.
11. **OutL**: Length of Outgoing part. If no outgoing part NaN is given.
12. **TotS**: Average step of full path.
13. **TotStdS**: Standard deviation of **TotS**.
14. **InpS**: Average step of Incoming part. If no incoming part NaN is given.
15. **InpStdS**: Standard deviation of **InpS**.
16. **ObjS**: Average step of Object part.
17. **ObjStdS**: Standard deviation of **ObjS**.
18. **OutS**: Average step of Outgoing part. If no outgoing part NaN is given.
19. **OutStdS**: Standard deviation of **OutS**.
20. **CType**: Cluster type of separate path.


## Separate path ID

Separate Path IDs are composed of three numbers separated by colon. First number is the layer number, if no *sandwich* option is used it is set to 0. The second number is residue number. Third number is consecutive number of the separate path made by the residue. Numeration starts with 0.

### 2.2.6.5 Frames dependent analysis

In addition to general summary Aqua-Duct calculates frames dependent parameters. Two types of values are calculated: number of traced paths, and *Object* and *Scope* sizes. Results are saved in the additional CSV file or are printed on the screen.

Calculated numbers of traced paths can be used to visualize behavior of the system in question. For example, one can analyze number of paths in two different clusters:



The above plot shows number of water molecules (or paths) in cluster 16 and 6 throughout the simulation. One can observe that number of molecules in cluster 6 diminishes approximately in the middle. This kind of plot can be easily generated with additional CSV data.

## Number of traced paths

For each frame, numbers of traced paths are calculated for following categories:

1. Name of traced molecules - `amol` is used for all possible names.
2. Paths types (object for standard paths and passing for passing paths) - `apaths` is used for all possible paths types.
3. Clusters and cluster types - `aclusts` is used for all possible clusters and `actypes` is used for all possible cluster types.
4. Part of paths. Possible values are: `walk`, `in`, `object`, `out`, and `in_out`. Where `walk` corresponds to any part of path and in case of passing paths only this category is used; `in`, `object`, and `out` correspond to incoming, object, and outgoing parts; `in_out` corresponds to sum of incoming and outgoing parts.

All the above listed categories are combined together, and the final number of calculated categories may be quite big.

## Size of *Object* and *Scope*

If option `calculate_scope_object_size` is set `True` and values of `scope_chull` and `object_chull` correspond to appropriate molecular entities, Aqua-Duct calculates area and volume of *Scope*

and *Object*. Calculated sizes are estimates as resulting from convex hull approximations.

## 2.2.7 Visualization

Sixth stage of *Valve* calculations visualizes results calculated in stages 1 to 4. Visualization is done with PyMOL. *Valve* creates visualizations in two modes:

1. Two files are created: special Python script and archive with data. Python script can be simply started with python, it automatically opens PyMol and loads all data from the archive. Optionally it can automatically save PyMol session.
2. PyMol is automatically started and all data is loaded directly to PyMol workspace.

Molecule is loaded as PDB file. Other objects like Inlets clusters or paths are loaded as CGO objects.

### 2.2.7.1 Visualization script

By default *Valve* creates Python visualization script and archive with data files. This script is a regular Python script. It does not depend on AQUA-DUCT. To run it, python2.7 and PyMol is required. If no **save** option is used *Valve* saves visualization script as `6_visualize_results.py`. To load full visualization call:

```
python 6_visualize_results.py --help

usage: 6_visualize_results.py [-h] [--save-session SESSION]
                             [--discard DISCARD] [--keep KEEP]
                             [--force-color FC] [--fast]

Aqua-Duct visualization script

optional arguments:
  -h, --help            show this help message and exit
  --save-session SESSION
                        Pymol session file name.
  --discard DISCARD     Objects to discard.
  --keep KEEP           Objects to keep.
  --force-color FC      Force specific color.
  --fast               Disable all objects while loading.
```

Option `--save-session` allows to save PyMol session file. Once visualization is loaded session is saved and PyMol closes. Option `--fast` increases slightly loading of objects.

Option `--force-color` allows to change default color of objects. It accepts list of specifications comprised of pairs 'object name' and 'color name'. For example: 'scope\_shape0 yellow cluster\_1 blue'. This will color *scope\_shape0* object in yellow and *cluster\_1* in blue:

```
python 6_visualize_results.py --force-color 'scope_shape0 yellow cluster_1 blue'
```

---

**Note:** List of specifications has to be given in parentheses.

---



---

**Note:** List of specifications has to comprise of full objects' names.

---



---

**Note:** Currently, `--force-color` does not allow to change color of molecules. It can be done in PyMol.

---

Options `--keep` and `--discard` allows to select specific objects for visualization. Both accept list of names comprising of full or partial object names. Option `--keep` instructs script to load only specified objects, whereas,

--discard instructs to skip specific objects. For example to keep shapes of object and scope, molecule and clusters only one can call following:

```
python 6_visualize_results.py --keep 'shape molecule cluster'
```

To discard all raw paths:

```
python 6_visualize_results.py --discard 'raw'
```

Options can be used simultaneously, order does matter:

1. If --keep is used first, objects are not displayed if they are not on the *keep* list. If they are on the list, visualization script checks if they are on the *discard* list. If yes, objects are not displayed.
2. If --discard is used first, objects are not displayed if they are on the *discard* list and are not on the *keep* list.

For example, in order to display molecule, clusters, and only raw master paths, one can use following command:

```
python 6_visualize_results.py --keep 'molecule cluster master' --discard 'smooth'
```

---

**Note:** Options --keep and --discard accepts both full and partial object names.

---

---

**Note:** List of names has to be given in parentheses.

---

### 2.2.7.2 Visualization objects

Following is a list of objects created in PyMOL (all of them are optional). PyMOL object names given in **bold** text or short explanation is given.

- Selected frame of the simulated system. Object name: **molecule** plus number of layer, if no *sandwich* option is used it becomes, by default, **molecule0**.
- Approximate shapes of object and scope. Objects names **object\_shape** and **scope\_shape** plus number of layer, if no *sandwich* option is used **0** is added by default.
- Inlets clusters, each cluster is a separate object. Object name: **cluster\_** followed by cluster annotation: outliers are annotated as **out**; regular clusters by ID.
- List of cluster types, raw paths. Each cluster type is a separate object. Object name composed of cluster type (colon replaced by underline) plus **\_raw**.
- List of cluster types, smooth paths. Each cluster type is a separate object. Object name composed of cluster type (colon replaced by underline) plus **\_smooth**.
- All raw paths. They can be displayed as one object or separated in to Incoming, Object and Outgoing part. Object name: **all\_raw**, or **all\_raw\_in**, **all\_raw\_obj**, and **all\_raw\_out**.
- All raw paths inlets arrows. Object name: **all\_raw\_paths\_io**.
- All smooth paths. They can be displayed as one object or separated in to Incoming, Object and Outgoing part. Object name: **all\_smooth**, or **all\_smooth\_in**, **all\_smooth\_obj**, and **all\_smooth\_out**.
- All raw paths inlets arrows. Object name: **all\_raw\_paths\_io**.
- Raw paths displayed as separate objects or as one object with several states. Object name: **raw\_paths\_** plus number of path or **raw\_paths** if displayed as one object.
- Smooth paths displayed as separate objects or as one object with several states. Object name: **smooth\_paths\_** plus number of path or **smooth\_paths** if displayed as one object.



- Raw paths arrows displayed as separate objects or as one object with several states. Object name: **raw\_paths\_io\_** plus number of path or **raw\_paths\_io** if displayed as one object.
- Smooth paths arrows displayed as separate objects or as one object with several states. Object name: **smooth\_paths\_io\_** plus number of path or **smooth\_paths\_io** if displayed as one object.

### 2.2.7.3 Color schemes

Inlets clusters are colored automatically. Outliers are gray.

Incoming parts of paths are red, Outgoing parts are blue. Object parts in case of smooth paths are green and in case of raw paths are green if residue is precisely in the object area or yellow if it leaved object area but it is not in the Outgoing part yet. *Passing paths* are displayed in grey.

Arrows are colored in accordance to the colors of paths.



## CONFIGURATION FILE OPTIONS

Valve configuration file is a simple and plain text file. It has similar structure as INI files commonly used in one of the popular operating systems and is compliant with Python module `ConfigParser`.

Configuration file comprises of several *sections*. They can be grouped into three categories. Names of sections are in **bold** text.

1. **Global settings:**

- **global**

2. **Stages options:**

1. **traceable\_residues**
2. **raw\_paths**
3. **separate\_paths**
4. **inlets\_clusterization**
5. **analysis**
6. **visualize**

3. **Methods options:**

- **smooth**
- **clusterization**
- **reclusteriation**

### 3.1 Section global

This section allows settings of trajectory data and is reserved for other future global options.

Option	Default value	Description
top	None	Path to topology file. Aqua-Duct supports PDB, PRMTOP, PFS topology files.
trj	None	Path to trajectory file. Aqua-Duct supports NC and DCD trajectory files.

Option **trj** can be used to provide list of trajectory files separated by standard path separator ':' on POSIX platforms and ';' on Windows - see `os.pathsep`.

---

**Note:** Options **top** and **trj** are mandatory.

---

## 3.2 Common settings of stage sections

Stages 1-4 which perform calculations have some common options allowing for execution control and saving/loading data.

Option	Default value	Description
execute	runonce	Option controls stage execution. It can have one of three possible values: <code>run</code> , <code>runonce</code> , and <code>skip</code> . If it is set to <code>run</code> calculations are always performed and if <b>dump</b> is set dump file is saved. If it is set to <code>runonce</code> calculations are performed if there is no dump file specified by <b>dump</b> option. If it is present calculations are skipped and data is loaded from the file. If it is set to <code>skip</code> calculations are skip and if <b>dump</b> is set data is loaded from the file.
dump	[dump file name]	File name of dump data. It is used to save results of calculations or to load previously calculated data - this depends on <b>execute</b> option. Default value of this option depends on the stage and for stages 1 to 4 is one of the following (listed in order): <ul style="list-style-type: none"> <li>• 1_traceable_residues_data.dump</li> <li>• 2_raw_paths_data.dump</li> <li>• 3_separate_paths_data.dump</li> <li>• 4_inlets_clusterization_data.dump</li> </ul>

Stages 5-6 also uses **execute** option, however, since they do not perform calculations *per se* in stead of **dump** option they use **save**.

Option	Default value	Description
execute	run	Option controls stage execution. It can have one of three possible values: <code>run</code> , <code>runonce</code> , and <code>skip</code> . If it is set to <code>run</code> or <code>runonce</code> stage is executed and results is saved according to <b>save</b> option. If it is set to <code>skip</code> stage is skipped.
save	[save file name]	File name for saving results. Default value of this option depends on the stage and for stages 1 to 4 is one of the following (listed in order): <ul style="list-style-type: none"> <li>• 5_analysis_results.txt &amp; 5_analysis_results.txt.csv</li> <li>• 6_visualize_results.py &amp; 6_visualize_results.tar.gz</li> </ul> Stage 5 saves <code>.txt</code> file with analysis results and, if requested, it saves additional <code>.csv</code> with various counts of traced molecules. Stage 6 can save results in two different ways: <ol style="list-style-type: none"> <li>1. As Python script - extension <code>.py</code> plus companion archive <code>.tar.gz</code>,</li> <li>2. As PyMOL session - extension <code>.pse</code>.</li> </ol>

## 3.3 Stage traceable\_residues

Option	Default value	Description
scope	None	Definition of <i>Scope</i> of interest. See also <a href="#">Scope definition</a> .
scope_convexhull	True	Flag to set if <i>Scope</i> is direct or convex hull definition.
scope_everyframe	False	Flag to set <i>Scope</i> evaluation mode. If set <code>True</code> <i>Scope</i> is evaluated in every frame. This make sense if the definition is complex and depends on distances between molecular entities.
object	None	Definition of <i>Object</i> of interest. See also <a href="#">Object definition</a> .

**Note:** Options **scope** and **object** are mandatory.

### 3.4 Stage raw\_paths

This stage also requires definition of the *Scope* and *Object*. If appropriate settings are not given, settings from the previous stage are used.

Option	Default value	Description
scope	None	Definition of <i>Scope</i> of interest. See also <i>Scope definition</i> . If <code>None</code> value from previous stage is used.
scope_convexhull	None	Flag to set if the <i>Scope</i> is direct or convex hull definition.
scope_everyframe	False	Flag to set <i>Scope</i> evaluation mode. If set <code>True</code> <i>Scope</i> is evaluated in every frame. This make sense if the definition is complex and depends on distances between molecular entities. If <code>None</code> value from previous stage is used.
object	None	Definition of <i>Object</i> of interest. See also <i>Object definition</i> . If <code>None</code> value from the previous stage is used
clear_in_object_info	False	If it is set to <code>True</code> information on occupation of <i>Object</i> site by traceable residues calculated in the previous stage is cleared and have to be recalculated. This is useful if definition of <i>Object</i> was changed.

### 3.5 Stage `separate_paths`

Option	Default value	Description
<code>discard_empty_paths</code>	True	If set to <code>True</code> empty paths are discarded.
<code>sort_by_id</code>	True	If set to <code>True</code> separate paths are sorted by ID. Otherwise they are sorted in order of appearance.
<code>discard_short_paths</code>	20	This option allows to discard paths which are shorter than the threshold which is defined as total number of frames.
<code>discard_short_object</code>	2.0	This option allows to discard paths which objects are shorter than the threshold which is defined as total length in metric units.
<code>discard_short_logic</code>	or	If both <code>discard_short_paths</code> and <code>discard_short_object</code> options are used, this option allows to set combination logic. If it is set <code>or</code> a path is discarded if any of discard criterion is met. If it is set <code>and</code> both criteria have to be met to discard path.
<code>auto_barber</code>	None	This option allows to select molecular entity used in Auto Barber procedure. See also <i>Auto Barber</i> and <i>barber_with_spheres()</i> .
<code>auto_barber_mincut</code>	None	Minimal radius of spheres used in Auto Barber. If a sphere has radius smaller then this value it is not used in AutoBarber procedure. This option can be switched off by setting it to <code>None</code> .
<code>auto_barber_maxcut</code>	2.8	Maximal radius of spheres used in Auto Barber. If a sphere has radius greater then this value it is not used in AutoBarber procedure. This option can be switched off by setting it to <code>None</code> .
<code>auto_barber_mincut_level</code>	True	If set <code>True</code> spheres of radius smaller than <b>mincut</b> are resized to <b>mincut</b> value.
<code>auto_barber_maxcut_level</code>	True	If set <code>True</code> spheres of radius greater than <b>maxcut</b> are resized to <b>maxcut</b> value.
<code>auto_barber_tovdw</code>	True	Correct cutting sphere by decreasing its radius by VdW radius of the closest atom.
<code>allow_passing_paths</code>	False	If set <code>True</code> paths that do not enter the object are detected and added to the rest of paths as 'passing' paths.

### 3.6 Stage inlets\_clusterization

Option	Default value	Description
recluster_outliers	False	If set to <code>True</code> reclusterization of outliers is executed according to the method defined in <b>reclusterization</b> section.
detect_outliers	False	If set, detection of outliers is executed. It could be set as a floating point distance threshold or set to <code>Auto</code> . See <i>Clusterization of inlets</i> for more details.
singletons_outliers	False	Maximal size of cluster to be considered as outliers. If set to number > 0 clusters of that size are removed and their objects are moved to outliers. See <i>Clusterization of inlets</i> for more details.
max_level	5	Maximal number of recursive clusterization levels.
create_master_paths	False	If set to <code>True</code> master paths are created (fast CPU and big RAM recommended; 50k frames long simulation may need ca 20GB of memory)
exclude_passing_in_clusterization	True	If set to <code>True</code> passing paths are not clustered with normal paths.
add_passing_to_clusters	None	Allows to run procedure for adding passing paths inlets to clusters with Auto Barber method. To enable this the option should be set to molecular entity that will be used by Auto Barber.

### 3.7 Stage analysis

Option	Default value	Description
dump_config	True	If set to <code>True</code> configuration options, as seen by Valve, are added to the head of results.
calculate_scope_object_size	False	If set to <code>True</code> volumes and areas of object and scope approximated by convex hulls will be calculated for each analyzed frames and saved in output CSV file.
scope_chull	None	Scope convex hull definition used in calculating volume and area.
object_chull	None	Object convex hull definition used in calculating volume and area.

### 3.8 Stage visualize

Option	Default value	Description
all_paths_raw	False	If <code>True</code> produces one object in PyMOL that holds all paths visualized by raw coordinates.
all_paths_smooth	False	If <code>True</code> produces one object in PyMOL that holds all paths visualized by smooth coordinates.
all_paths_split	False	If is set <code>True</code> objects produced by <b>all_paths_raw</b> and <b>all_paths_smooth</b> are split into Incoming, Object, and Outgoing parts and visualized as three different objects.
all_paths_raw_io	False	If set <code>True</code> arrows pointing beginning and end of paths are displayed oriented accordingly to raw paths orientation.

Continued on next page

Table 1 – continued from previous page

Option	Default value	Description
<code>all_paths_smooth_io</code>	False	If set True arrows pointing beginning and end of paths are displayed oriented accordingly to smooth paths orientation.
<code>simply_smooths</code>	RecursiveVector	<p>Option indicates linear simplification method to be used in plotting smooth paths. Simplification removes points which do not (or almost do not) change the shape of smooth path. Possible choices are:</p> <ul style="list-style-type: none"> <li>RecursiveVector (<code>LinearizeRecursiveVector</code>),</li> <li>HobbitVector (<code>LinearizeHobbitVector</code>),</li> <li>OneWayVector (<code>LinearizeOneWayVector</code>),</li> <li>RecursiveTriangle (<code>LinearizeRecursiveTriangle</code>),</li> <li>HobbitTriangle (<code>LinearizeHobbitTriangle</code>),</li> <li>OneWayTriangle (<code>LinearizeOneWayTriangle</code>).</li> </ul> <p>Optionally name of the method can be followed by a threshold value in parentheses, i.e. <code>RecursiveVector(0.05)</code>. For sane values of thresholds see appropriate documentation of each method. Default values work well. This option is not case sensitive. It is recommended to use default method or <code>HobbitVector</code> method.</p>
<code>paths_raw</code>	False	If set True raw paths are displayed as separate objects or as one object with states corresponding to number of path.
<code>paths_smooth</code>	False	If set True smooth paths are displayed as separate objects or as one object with states corresponding to number of path.
<code>paths_raw_io</code>	False	If set True arrows indicating beginning and end of paths, oriented accordingly to raw paths, are displayed as separate objects or as one object with states corresponding to number of paths.
<code>paths_smooth_io</code>	False	If set True arrows indicating beginning and end of paths, oriented accordingly to smooth paths, are displayed as separate objects or as one object with states corresponding to number of paths.
<code>paths_states</code>	False	If True objects displayed by <b><code>paths_raw</code></b> , <b><code>paths_smooth</code></b> , <b><code>paths_raw_io</code></b> , and <b><code>paths_smooth_io</code></b> are displayed as one object with states corresponding to number of paths. Otherwise they are displayed as separate objects.
<code>ctypes_raw</code>	False	Displays raw paths in a similar manner as non split <b><code>all_paths_raw</code></b> but each cluster type is displayed in separate object.
<code>ctypes_smooth</code>	False	Displays smooth paths in a similar manner as non split <b><code>all_paths_smooth</code></b> but each cluster type is displayed in separate object.
<code>show_molecule</code>	False	If is set to selection of some molecular object in the system, for example to <code>protein</code> , this object is displayed.
<code>show_molecule_frames</code>	0	Allows to indicate which frames of object defined by <b><code>show_molecule</code></b> should be displayed. It is possible to set several frames. In that case frames would be displayed as states.

Continued on next page



Table 1 – continued from previous page

Option	Default value	Description
<code>show_scope_chull</code>	False	If is set to selection of some molecular object in the system, for example to <code>protein</code> , convex hull of this object is displayed.
<code>show_scope_chull_frames</code>	0	Allows to indicate for which frames of object defined by <b>show_chull</b> convex hull should be displayed. It is possible to set several frames. In that case frames would be displayed as states.
<code>show_object_chull</code>	False	If is set to selection of some molecular object in the system convex hull of this object is displayed. This works exactly the same way as <b>show_chull</b> but is meant to mark object shape. It can be achieved by using <i>name</i> * <i>and</i> molecular object definition plus some spatial constraints, for example those used in object definition.
<code>show_object_chull_frames</code>	0	Allows to indicate for which frames of object defined by <b>show_object</b> convex hull should be displayed. It is possible to set several frames. In that case frames would be displayed as states.

---

**Note:** Possibly due to limitations of MDAnalysis only whole molecules can be displayed. If **show\_molecule** is set to `backbone` complete protein will be displayed any way. This may change in future version of MDAnalysis and or *aqueduct*.

---



---

**Note:** If several frames are selected they are displayed as states which may interfere with other PyMOL objects displayed with several states.

---



---

**Note:** If several states are displayed protein tertiary structure data might be lost. This seems to be limitation of either MDAnalysis or PyMOL.

---

## 3.9 Clusterization sections

Default section for definition of clusterization method is named **clusterization** and default section for reclusterization method definition is named **reclusterization**. All clusterization sections shares some common options. Other options depends on the method.

Option	Default value	Description
method	barber or dbscan	Name of clusterization method. It has to be one of the following: barber, dbscan, affprop, meanshift, birch, kmeans. Default value depends whether it is <b>clusterization</b> section (barber) or <b>reclusterization</b> section (dbscan).
recursive_clusterization	clusterization or None	If it is set to name of some section that holds clusterization method settings this method will be called in the next recursion of clusterization. Default value for <b>reclusterization</b> is None.
recursive_threshold	None	Allows to set threshold that excludes clusters of certain size from reclustering. Value of this option comprises of <i>operator</i> and <i>value</i> . Operator can be one of the following: >, >=, <=, <. Value have to be expressed as floating number and it have to be in the range of 0 to 1. One can use several definitions separated by a space character. Only clusters of size complying with all thresholds definitions are submitted to reclustering.

### 3.9.1 barber

Clusterization by **barber** method bases on *Auto Barber* procedure. For each inlets a sphere is constructed according to Auto Barber **separate\_paths** stage settings or according to parameters given in clusterization section. Next, inlets that form coherent clouds of mutually intersecting spheres are grouped in to clusters. Method **barber** supports the same settings as Auto Barber settings:

Option	Value type	Description
auto_barber	str	This option allows to select molecular entity used in Auto Barber procedure. See also <i>Auto Barber</i> and <i>barber_with_spheres()</i> .
auto_barber_mincut	float	Minimal radius of spheres used in Auto Barber. If a sphere has radius smaller then this value it is not used to cut. This option can be switched off by setting it to <i>None</i> .
auto_barber_maxcut	float	Maximal radius of spheres used in Auto Barber. If a sphere has radius greater then this value it is not used to cut. This option can be switched off by setting it to <i>None</i> .
auto_barber_mincut_level	bool	If set <i>True</i> spheres of radius less then <b>mincut</b> are resized to <b>mincut</b> value.
auto_barber_maxcut_level	bool	If set <i>True</i> spheres of radius greater then <b>maxcut</b> are resized to <b>maxcut</b> value.
auto_barber_tovdw	bool	Correct cutting sphere by decreasing its radius by VdW radius of the closest atom.

### 3.9.2 dbscan

For detailed description look at `sklearn.cluster.DBSCAN` documentation. Following table summarized options available in *Valve* and is a copy of original documentation.

Option	Value type	Description
eps	float	The maximum distance between two samples for them to be considered as in the same neighborhood.
min_samples	int	The number of samples (or total weight) in a neighborhood for a point to be considered as a core point. This includes the point itself.
metric	str	The metric to use when calculating distance between instances in a feature array. Can be one of the following: <ul style="list-style-type: none"> <li>• euclidean,</li> <li>• cityblock,</li> <li>• cosine,</li> <li>• manhattan.</li> </ul>
algorithm	str	The algorithm to be used by the NearestNeighbors module to compute pointwise distances and find nearest neighbors. Can be one of the following: <ul style="list-style-type: none"> <li>• auto,</li> <li>• ball_tree,</li> <li>• kd_tree,</li> <li>• brute.</li> </ul>
leaf_size	int	Leaf size passed to BallTree or cKDTree.

### 3.9.3 affprop

For detailed description look at [AffinityPropagation](#) documentation. Following table summarized options available in *Valve* and is a copy of original documentation.

Option	Value type	Description
damping	float	Damping factor between 0.5 and 1.
convergence_iter	int	Number of iterations with no change in the number of estimated clusters that stops the convergence.
max_iter	int	Maximum number of iterations.
preference	float	Points with larger values of preferences are more likely to be chosen as exemplars.

### 3.9.4 meanshift

For detailed description look at [MeanShift](#) documentation. Following table summarized options available in *Valve* and is a copy of original documentation.

Option	Value type	Description
bandwidth	Auto or float	Bandwidth used in the RBF kernel. If Auto or None automatic method for bandwidth estimation is used. See <a href="#">estimate_bandwidth()</a> .
cluster_all	bool	If true, then all points are clustered, even those orphans that are not within any kernel.
bin_seeding	bool	If true, initial kernel locations are not locations of all points, but rather the location of the discretized version of points, where points are binned onto a grid whose coarseness corresponds to the bandwidth.
min_bin_freq	int	To speed up the algorithm, accept only those bins with at least min_bin_freq points as seeds. If not defined, set to 1.

### 3.9.5 birch

For detailed description look at [Birch](#) documentation. Following table summarized options available in *Valve* and is a copy of original documentation.

Option	Value type	Description
threshold	float	The radius of the subcluster obtained by merging a new sample and the closest subcluster should be lesser than the threshold. Otherwise a new subcluster is started.
branching_factor	int	Maximum number of CF subclusters in each node.
n_clusters	int	Number of clusters after the final clustering step, which treats the subclusters from the leaves as new samples. By default, this final clustering step is not performed and the subclusters are returned as they are.

### 3.9.6 kmeans

For detailed description look at [KMeans](#) documentation. Following table summarized options available in *Valve* and is a copy of original documentation.

Option	Value type	Description
n_clusters	int	The number of clusters to form as well as the number of centroids to generate.
max_iter	int	Maximum number of iterations of the k-means algorithm for a single run.
n_init	int	Number of time the k-means algorithm will be run with different centroid seeds. The final results will be the best output of n_init consecutive runs in terms of inertia.
init	str	Method for initialization, defaults to k-means++. Can be one of following: k-means++ or random.
tol	float	Relative tolerance with regards to inertia to declare convergence.

## 3.10 Smooth section

Section **smooth** supports following options:

Option	Value type	Description
method	str	Smoothing method. Can be one of the following: <ul style="list-style-type: none"><li>• window, (see <a href="#">WindowSmooth</a>)</li><li>• mss, (see <a href="#">MaxStepSmooth</a>)</li><li>• window_mss, (see <a href="#">WindowOverMaxStepSmooth</a>)</li><li>• awin, (see <a href="#">ActiveWindowSmooth</a>)</li><li>• awin_mss, (see <a href="#">ActiveWindowOverMaxStepSmooth</a>)</li><li>• dwin, (see <a href="#">DistanceWindowSmooth</a>)</li><li>• dwin_mss, (see <a href="#">DistanceWindowOverMaxStepSmooth</a>)</li><li>• savgol. (see <a href="#">SavgolSmooth</a>)</li></ul>
recursive	int	Number of recursive runs of smoothing method.
window	int or float	In window based method defines window size. In plain window it has to be int number. In savgol it has to be odd integer.
step	int	In step based method defines size of the step.
function	str	In window based methods defines averaging function. Can be mean or median.
polyorder	int	In savgol is polynomial order.

## VALVE TUTORIAL

This tutorial assumes *aqueduct* and *Valve* is already installed - see *Aqua-Duct installation guide*. It is also assumed that user is acquainted with *Valve manual* and *Valve Configuration file options*.

### 4.1 Valve invocation

Usually *Valve* is run by:

```
valve.py
```

To check if *Valve* is installed and works properly try to issue following commands:

```
valve.py --help  
valve.py --version
```

### 4.2 Test data

#### Mouse!

We will use 1 ns MD simulation data of sEH protein (PDBID **1cqz**). This simulation was performed in Amber 14. Necessary files can be found at [Aqua-Duct home page](#) in section [download](#). Required data is in the *sample data* file.

### 4.3 Inspect your system

Before we start any calculations let's have a look at the protein of interest. Start PyMOL and get 1cqz PDB structure (for example by typing in PyMOL command prompt `fetch 1cqz`).

To setup *Valve* calculations we need to know the active site of the protein. More precisely we need to know IDs of residues that are in the active site. This would allow us to create *Object definition*.

But wait. Is it really the correct structure? How many chains there are? What is the numeration of residues? How does it compare with the topology file from *sample data*?

#### 4.3.1 Create *Object definition*

Let's load another structure. Open file `1cqz_sample_topology.pdb` (see *Test data*). It is a first frame of the MD simulation and it is an example of how the frame of MD looks like. In order to create *Object definition* you have to discover following things:

1. What is the name of water molecules?
2. What are numbers of residues in the active site?

3. What size the active site is of?

---

**Note:** It is also a good idea to open `.pdb` file in your favorite text editor and look at residue numbers and names.

---

### 4.3.2 Create *Scope definition*

*Scope definition* is easy to create. We will use *Convex hull* version so the scope definition could be simply backbone.

## 4.4 Prepare config file

*Valve* performs calculations according to the configuration (aka *config*) file.

Lets start from dumping config file template to `config.txt` file. Open it in your favorite editor and fill all options. If you have troubles look at *Configuration file options* (and *Valve manual*).

Things to remember:

1. Provide correct paths to topology and trajectory data.
2. Enter correct *Object* and *Scope* definitions.
3. Make sure visualization is switched on.

## 4.5 Run *Valve*

Make sure all necessary data is in place. Open terminal, go to your working directory and type in:

```
valve.py -c config.txt
```

Depending on your machine and current load it may take a while (matter of minutes) to complete all calculations.

### 4.5.1 Visual inspection

In the last stage PyMOL should pop up and *Valve* should start to feed it with visualization data. This would take a moment and if you set up `save` option a PyMOL session would be saved. Once it is done *Valve* quits and switches off PyMOL. Now, you can restart it and read saved session.

### 4.5.2 Clusterization

Improve clusterization of Inlets. See *Configuration file options* for more hints on available clusterization options.

### 4.5.3 Analysis tables

Open `5_analysis_results.txt` file and look at summaries and tables. See also *Valve manual*.

## 4.6 Feedback

Give us your opinion. Send your questions, inquiries, anything to developer(s): [info@aquaduct.pl](mailto:info@aquaduct.pl). There are couple of questions that might be useful to form your opinion.

1. What do you like in *Valve* and *Aqua-Duct*?
2. What do you do not like in *Valve* or *Aqua-Duct*?
3. What is missing?
4. Do you find it useful?





## AQUADUCT

## 5.1 aquaduct package

### 5.1.1 Subpackages

#### 5.1.1.1 aquaduct.apps package

##### Submodules

##### aquaduct.apps.data module

**class** GlobalConfigStore

Bases: `object`

**cachedir** = None

**cachemem** = False

**class** CoordsRangeIndexCache

Bases: `object`

**cache** = {}

**get\_cric\_reader** (*mode*='r')

**save\_cric** ()

**load\_cric** ()

**check\_version\_compliance** (*current*, *loaded*, *what*)

**check\_versions** (*version\_dict*)

**class** LoadDumpWrapper (*filehandle*)

Bases: `object`

This is wrapper for pickled data that provides compatibility with earlier versions of Aqua-Duct.

Conversions in use:

- 1) replace 'aquaduct.' by 'aquaduct.'

**\_\_init\_\_** (*filehandle*)

x.**\_\_init\_\_**(...) initializes x; see help(type(x)) for signature

**convert** (*s*)

**read** (*\*args*, *\*\*kwargs*)

**readline** (*\*args*, *\*\*kwargs*)

**get\_vda\_reader** (*filename*)

```
class ValveDataAccess_pickle (mode=None, data_file_name=None, reader=None)
    Bases: aquaduct.apps.data.ValveDataAccess
    mimic_old_var_name = 'aq_data_to_save'
    unknown_names = 'UNK'
    open (data_file_name, mode)
    close ()
    load ()
    dump (**kwargs)
    get_variable (name)
    set_variable (name, value)

class ValveDataAccessRoots
    Bases: object
    roots = []
    open (data_file_name, mode)
    close_all ()
    __del__ ()

get_object_name (something)
get_object_from_name (name)

class IdsOverIds
    Bases: object
    static dict2arrays (d)
    static arrays2dict (values=None, keys_lens=None)

class ValveDataAccess_nc (*args, **kwargs)
    Bases: aquaduct.apps.data.ValveDataAccess
    __init__ (*args, **kwargs)
        x.__init__(...) initializes x; see help(type(x)) for signature
    open (data_file_name, mode)

ValveDataAccess
    alias of aquaduct.apps.data.ValveDataAccess_pickle
```

## **aquaduct.apps.valvecore module**

```
class ValveConfig
    Bases: object, aquaduct.apps.valvecore.ConfigSpecialNames
    __init__ ()
        x.__init__(...) initializes x; see help(type(x)) for signature
    static common_config_names ()
    static common_traj_data_config_names ()
    static global_name ()
    static cluster_name ()
    static recluster_name ()
    static recursive_clusterization_name ()
```

```

static recursive_threshold_name ()
static smooth_name ()
stage_names (nr=None)
get_common_traj_data (stage)
get_global_options ()
get_stage_options (stage)
get_cluster_options (section_name=None)
get_recluster_options ()
get_smooth_options ()
get_default_config ()
load_config (filename)
save_config_stream (fs)
save_config (filename)
get_general_comment (section)
dump_config (dump_template=False)
_ValveConfig__make_options_nt (input_options)
valve_begin ()
valve_end ()
valve_load_config (filename, config)
valve_exec_stage (stage, config, stage_run, no_io=False, run_status=None, **kwargs)
stage_I_run (config, options, **kwargs)
stage_II_run (config, options, all_res=None, number_frame_rid_in_object=None, **kwargs)
stage_III_run (config, options, paths=None, **kwargs)
stage_IV_run (config, options, spaths=None, center_of_system=None, **kwargs)
stage_V_run (config, options, spaths=None, paths=None, inls=None, ctypes=None, reader=None,
             **kwargs)
stage_VI_run (config, options, spaths=None, inls=None, ctypes=None, master_paths=None, mas-
              ter_paths_smooth=None, **kwargs)
aquaduct_version_nice ()
    Returns aquaduct version number as nicely formatted string.
        Returns string composed on the basis of the number returned by version().
        Return type str

```

## Module contents

### 5.1.1.2 aquaduct.geom package

#### Submodules

#### aquaduct.geom.cluster module

This module provides functions for clusterization. Clusterization is done by `scikit-learn` module.

**get\_required\_params** (*method*)

**class BarberClusterResult** (*labels\_*)

Bases: `object`

Helper class for results of barber clusterization.

**\_\_init\_\_** (*labels\_*)

`x.__init__(...)` initializes `x`; see `help(type(x))` for signature

**class BarberCluster**

Bases: `object`

Wrapper class that implements *barber* clusterization.

**fit** (*coords*, *spheres=None*)

**Parameters**

- **coords** (*Iterable*) – Input coordinates of points to be clustered.
- **spheres** (*Iterable*) – Input spheres for each point.

**MeanShiftBandwidth** (*X*, *\*\*kwargs*)

Helper function for automatic calculation of a bandwidth for MeanShift method.

**Parameters** **X** (*Iterable*) – Coordinates of points to be clustered.

**class PerformClustering** (*method*, *\*\*kwargs*)

Bases: `object`

Helper class for clusterization.

**\_\_init\_\_** (*method*, *\*\*kwargs*)

**Parameters** **method** (*object*) – Class that implements cclusterization via *fit* method.

**\_\_str\_\_** () <==> *str(x)*

**\_\_call\_\_** (...) <==> *x(...)*

**\_get\_noclusters** (*n*)

**fit** (*coords*, *spheres=None*)

**Parameters**

- **coords** (*Iterable*) – Input coordinates of points to be clustered.
- **spheres** (*Iterable*) – Input spheres for each point. Optional, important only if method is *BarberCluster*.

**Returns** Clusters numbers.

**Return type** list of int

**centers** ()

**Returns** Centers of clusters.

## **aquaduct.geom.convexhull module**

**\_vertices\_ids** (*convexhull*)

**\_vertices\_points** (*convexhull*)

**\_point\_within\_convexhull** (*convexhull*, *point*)

**\_facets** (*convexhull*)

**\_edges** (*\*args*, *\*\*kwargs*)

**is\_point\_within\_convexhull** (*point\_chull*)

**aquaduct.geom.master module**

**part2type\_dict** = {0: 's', 1: 'c', 2: 's'}

Part number to *GenericPathTypeCodes* dictionary.

**parts** = (0, 1, 2)

Parts enumerate.

**class CTypeSpathsCollectionWorker** (*spaths=None, ctype=None, bias\_long=5, smooth=None, lock=None*)

Bases: *object*

Worker class for averaging spaths in points of master path.

**\_\_init\_\_** (*spaths=None, ctype=None, bias\_long=5, smooth=None, lock=None*)

Core method for averaging spaths in to master path.

Averaging is done in chunks.

**Parameters**

- **spaths** (*list*) – List of separate paths to average.
- **ctype** (*InletClusterGenericType*) – CType of spaths.
- **bias\_long** (*int*) – Bias towards long paths used in *lens\_norm()*.
- **smooth** (*Smooth*) – Smoothing method.

**coords\_types\_prob\_widths** (*sp\_slices\_*)

Calculates average coordinates, type and width in given chunk.

Parameter *sp\_slices\_* is tuple of length equal to number of spaths. It contains slices for all spaths respectively. With these slices spaths are cut and **only** resulting chunks are used for calculations.

Therefore, this method average spaths in one point of master math. This point is defined by slices submitted as *sp\_slices\_* parameter.

Algorithm of averaging (within current chunks of spaths):

1. Coordinates for all spaths are collected.
2. Lengths of all spaths are collected (from cached variables) and kept as lists of lengths equal to chunks' sizes.

---

**Note:** Lengths of collected lengths of spaths are of the same size as coordinates

---

3. New coordinates are calculated as weighted average of collected coordinates with *numpy.average()*. As weights collected lengths are used.

---

**Note:** Function *numpy.average()* is called with flatten coordinates and lengths.

---

4. Width of average path is calculated as mean value of flatten coordinates mutual distances.
5. Type of average paths is calculated as probability (frequency) of *scope\_name*.

**Parameters** *sp\_slices* (*tuple*) – Slices that cut chunks from all paths.

**Return type** 3 element tuple

**Returns** coordinates, type (frequency), and width of averaged spaths in current point

**\_\_call\_\_** (*nr\_sp\_slices\_*)

Callable interface.

**Parameters** *nr\_sp\_slices* (*tuple*) – Two element tuple: *nr* and *sp\_slice*

```
class CTypeSpathsCollection (spaths=None,    ctype=None,    bias_long=5,    pbar=None,
                             threads=1)
```

Bases: `object`

Object for grouping separate paths that belong to the same CType. Method `get_master_path()` allows for calculation of average path.

```
parts = (0, 1, 2)
```

Enumeration of spath parts.

```
__init__ (spaths=None, ctype=None, bias_long=5, pbar=None, threads=1)
```

#### Parameters

- **spaths** (*list*) – List of separate paths.
- **ctype** (`InletClusterGenericType`) – CType of spaths.
- **bias\_long** (*int*) – Bias towards long paths used in `lens_norm()`.
- **pbar** – Progress bar object.
- **threads** (*int*) – Number of available threads.

```
beat ()
```

Touch progress bar, if any.

```
update ()
```

Update progress bar by one, if any.

```
lens ()
```

Returns total lengths of all paths.

If ctype in `#:#` and not 0 and not None then take length of *object* part only.

**Returns** Total (or *object* part) lengths of all paths.

**Return type** `numpy.ndarray`

```
lens_norm ()
```

Returns normalized lengths calculated by `lens()`.

Applied normalization is twofold:

1. All lengths are divided by maximal length, and
2. All lengths are subjected to `pow()` function with `p = bias_long`.

**Returns** Normalized total (or *object* part) lengths of all paths.

**Return type** `numpy.ndarray`

```
lens_real ()
```

Returns real lengths of all paths.

**Returns** Sizes of all paths.

**Return type** `list`

```
full_size ()
```

Returns desired size of master path.

**Returns** Size of master path.

**Return type** `int`

```
static simple_types_distribution (types)
```

Calculates normalized sizes of incoming, object, and outgoing parts of spath using generic types.

It is assumed that spath has object part.

**Parameters** **types** (*list*) – List of generic types.

**Return type** 3 element list

**Returns** Normalized sizes of incomin, object, and outgoing parts.

**types\_distribution()**

**Return type** `numpy.matrix`

**Returns** median values of `simple_types_distribution()` for all spaths.

**types\_prob\_to\_types** (*types\_prob*)

Changes types probabilities as returned by `CTypeSpathsCollectionWorker.coords_types_prob_widths()` to types.

**Parameters** **types\_prob** (*list*) – List of types probabilities.

**Return type** list

**Returns** List of `GenericPathTypeCodes`.

**get\_master\_path** (*smooth=None, resid=(0, 0)*)

Averages spaths into one master path.

This is done in steps:

1. Master path is an average of bunch of spaths. Its length is determined by `full_size()` method.
2. All spaths are then divided in to chunks according to `xzip_xzip()` function with N set to lenght of master path. This results in list of length equal to the length of master path. Elements of this lists are slice objects that can be used to slice spaths in appropriate chunks.
3. Next, for each element of this list `CTypeSpathsCollectionWorker.coords_types_prob_widths()` method is called. Types probabilities are changed to types wiht `types_prob_to_types()`.
4. Finally, all data are used to create appropriate MasterPath. If this fails `None` is returned.

**Parameters**

- **smooth** (*Smooth*) – Smoothing method.
- **resid** (*int*) – Residue ID of master path.

**Return type** `MasterPath`

**Returns** Average path as `MasterPath` object or `None` if creation of master path failed.

**class FakeSingleResidueSelection** (*resid, frames, coords*)

Bases: `aquaduct.traj.sandwich.SingleResidueSelection`

**\_\_init\_\_** (*resid, frames, coords*)

`x.__init__(...)` initializes x; see `help(type(x))` for signature

**coords** (*\*\*kwargs*)

**coords\_smooth** (*sranges, smooth*)

## aquaduct.geom.pca module

**class Center** (*X*)

Bases: `object`

**\_\_init\_\_** (*X*)

`x.__init__(...)` initializes x; see `help(type(x))` for signature

**\_\_call\_\_** (*...*)  $\Leftrightarrow x(...)$

**undo** (*X*)

```
class Normalize (X)
    Bases: object

    __init__ (X)
        x.__init__(...) initializes x; see help(type(x)) for signature

    __call__ (...) <==> x(...)

    undo (X)

class Standartize (X)
    Bases: aquaduct.geom.pca.Center, aquaduct.geom.pca.Normalize

    __init__ (X)
        x.__init__(...) initializes x; see help(type(x)) for signature

    __call__ (...) <==> x(...)

    undo (X)

class PCA (X, prepro=None)
    Bases: object

    __init__ (X, prepro=None)
        x.__init__(...) initializes x; see help(type(x)) for signature

    P

    preprocess (X)

    preprocess_undo (X)

    __call__ (...) <==> x(...)

    undo (T)
```

## aquaduct.geom.smooth module

Smooth module defines methods for smoothing of trajectories.

Available methods:

<code>SavgolSmooth</code>	Savitzky-Golay based smoothing.
<code>WindowSmooth</code>	Defined size window smoothing.
<code>DistanceWindowSmooth</code>	Distance defined size window smoothing.
<code>ActiveWindowSmooth</code>	Active size window smoothing.
<code>MaxStepSmooth</code>	Maximal step smoothing.
<code>WindowOverMaxStepSmooth</code>	Window smoothing over maximal step smoothing.
<code>DistanceWindowOverMaxStepSmooth</code>	Distance window smoothing over maximal step smoothing.
<code>ActiveWindowOverMaxStepSmooth</code>	Active window smoothing over maximal step smoothing.

```
class Smooth (recursive=None, **kwargs)
    Bases: object

    Base class for all smoothing methods.

    __init__ (recursive=None, **kwargs)

        Parameters recursive (int) – Number of recursions of the method, everything evaluated to False is equivalent to 1.

    smooth (coords)
        Abstract method for smoothing method implementation.
```



**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

`__call__` (*coords*)

Call method for all smoothing methods.

Input coordinates should be iterable and each element should be `numpy.ndarray`. If length of `coords` is less then 3 smoothing method is not run and coordinates are returned unchanged.

If `recursive` is set smoothing method is applied appropriate number of times.

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**Return type** `numpy.ndarray`

**Returns** Smoothed coordinates.

**class** `GeneralWindow` (*function=<function mean>, \*\*kwargs*)

Bases: `object`

Base class for window based smoothing methods.

`__init__` (*function=<function mean>, \*\*kwargs*)

**Parameters** `function` (*function*) – Function to be used for averaging coordinates within a window.

**static** `max_window_at_pos` (*pos, size*)

Method returns maximal possible window at given position of the list with given size of the list. Returned window fits in to the list of given size and is symmetrical.

**Parameters**

- `pos` (*int*) – Position in question.
- `size` (*int*) – Length of the list.

**Return type** 2 element tuple of int

**Returns** Lowest possible bound and highest possible bound of the window.

**check\_bounds\_at\_max\_window\_at\_pos** (*lb, ub, pos, size*)

Method checks if window fits in to maximal possible window calculated according to `max_window_at_pos()`. If not window is corrected.

**Parameters**

- `lb` (*int*) – Lower bound of the window in question.
- `ub` (*int*) – Upper bound of the window in question.
- `pos` (*int*) – Position in question.
- `size` (*int*) – Length of the list.

**Return type** 2 element tuple of int

**Returns** Lowest possible bound and highest possible bound of the window corrected to maximal possible window.

**class** `IntWindow` (*window=5, \*\*kwargs*)

Bases: `aquaduct.geom.smooth.GeneralWindow`

Base class for all window smoothing methods that require integer window.

`__init__` (*window=5, \*\*kwargs*)

**Parameters** `window` (*int*) – One side size of the window.

**class** `FloatWindow` (*window=5.0, \*\*kwargs*)

Bases: `aquaduct.geom.smooth.GeneralWindow`

Base class for all window smoothing methods that require float window.

`__init__` (*window=5.0, \*\*kwargs*)

**Parameters** `window` (*float*) – Size of the window.

**class** `WindowSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`, `aquaduct.geom.smooth.IntWindow`

Defined size window smoothing.

For each coordinate a symmetrical (if possible) window of size defined by `window` is created. In case of coordinates at the edges created window is truncated to the edges. Next, all coordinates within the window are averaged with a function defined by `function`. Resulting value(s) are the smoothed coordinates.

`__init__` (*\*\*kwargs*)

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth` (*\*\*kwargs*)

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `DistanceWindowSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`, `aquaduct.geom.smooth.FloatWindow`

Distance defined size window smoothing.

This is modification of `WindowSmooth` method. The difference is in the definition of the window size. Here, it is an average distance between points of input coordinates. Thus, before smoothing average distance between all points is calculated and this value is used to calculate actual window size.

Next, for each coordinate a symmetrical (if possible) window of size calculated in the first step is created. In case of coordinates at the edges created window is truncated to the edges. Next, all coordinates within the window are averaged with a function defined by `function`. Resulting value(s) are the smoothed coordinates.

`__init__` (*\*\*kwargs*)

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth` (*\*\*kwargs*)

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `ActiveWindowSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`, `aquaduct.geom.smooth.FloatWindow`

Active size window smoothing.

Similarly to `DistanceWindowSmooth` method the window size is defined as a distance. The difference is that the actual window size is calculated for each point separately. Thus, for each coordinate the window is calculated by examining the distance differences between points. In this method window is not necessarily symmetrical. Once window is calculated all coordinates within the window are averaged with a function defined by `function`. Resulting value(s) are the smoothed coordinates.

`__init__` (*\*\*kwargs*)

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth` (*\*\*kwargs*)

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `MaxStepSmooth` (*step=1.0*, *\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`

Maximal step smoothing.

This method moves thorough coordinates and calculates distance over the traversed path. If it is then `step` the coordinate is used as a “cardinal point”. The beginning and the end of the path are also added to the list of cardinal points. Next, all cardinal points and points of linear interpolation between cardinal points

are returned as smoothed coordinates. Number of interpolated points is in accordance to points skipped between cardinal points.

`__init__ (step=1.0, **kwargs)`

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth (**kwargs)`

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `SavgolSmooth` (*window\_length=5, polyorder=2, \*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`

Savitzky-Golay based smoothing.

Method uses 1D filter available in SciPy, see `savgol_filter()`. For each dimension filter is applied separately. Only `window_length` and `polyorder` attributes are used.

`__init__ (window_length=5, polyorder=2, **kwargs)`

**Param** `int window_length`: Size of the window, odd number.

**Param** `int polyorder`: Polynomial order.

`set_savgol_function ()`

`smooth (**kwargs)`

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `WindowOverMaxStepSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`

Window smoothing over maximal step smoothing.

First, `MaxStepSmooth` is applied, and then `WindowSmooth`.

`__init__ (**kwargs)`

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth (coords)`

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `ActiveWindowOverMaxStepSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`

Active window smoothing over maximal step smoothing.

First, `MaxStepSmooth` is applied, and then `ActiveWindowSmooth`.

`__init__ (**kwargs)`

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth (coords)`

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

**class** `DistanceWindowOverMaxStepSmooth` (*\*\*kwargs*)

Bases: `aquaduct.geom.smooth.Smooth`

Distance window smoothing over maximal step smoothing.

First, `MaxStepSmooth` is applied, and then `DistanceWindowSmooth`.

`__init__ (**kwargs)`

**Parameters** `recursive` (*int*) – Number of recursions of the method, everything evaluated to `False` is equivalent to 1.

`smooth` (*coords*)

**Parameters** `coords` (*Iterable*) – Input coordinates to be smoothed.

## **aquaduct.geom.traces module**

`diff` (*trace*)

This function calculates the distance between 2 given points.

**Parameters** `trace` – coordinates in numpy array object

**Returns** distance between points

`tracepoints` (*start, stop, nr*)

**Parameters**

- **start** – coordinates of the first point as a numpy array object
- **stop** – coordinates of the second point as a numpy array object
- **nr** – number of elements between the first and second point

**Returns** two-dimentional numpy array; number of dimentions depends on nr parameter

`midpoints` (*paths*)

The function returns a tuple of numpy arrays extended with mid point spanning last and first element(column) of these arrays.

**Parameters** `paths` – a tuple of 2-dimentional np.arrays that hold 3D coordinates; each element holds one trace, all elements are supposed to make one path divided in to sections

**Returns** paths elements with additional mid points as a generator object

`length_step_std` (*trace*)

This function calculates sum, mean and standard deviation from all segments of a trace.

**Parameters** `trace` – coordinates of points as numpy array

**Returns** a tuple with basics statistics of a trace

`derrivative` (*values*)

`vector_norm` (*V*)

**Parameters** `V` – a vector in a form of array-like object, tuple or a list

**Returns** normalized length of a vector

`triangle_angles` (*A, B, C*)

Parameters are coordinates of points which are tops of triangle. The function calculates angles in a triangle formed by given coordinates.

**Parameters**

- **A** – coordinates of the first point
- **B** – coordinates of the second point
- **C** – coordinates of the third point

**Returns** list of arguments where angle is given in radians , the output is as follow:  
[BAC,CAB,ABC]

`triangle_angles_last` (*A, B, C*)

Parameters are coordinates of points which are tops of triangle. The function calculates the [ABC] angle.

**Parameters**

- **A** – coordinates of the first point [A top]
- **B** – coordinates of the second point [B top]
- **C** – coordinates of the third point [C top]

**Returns** list with one value of ABC angle in radians

**triangle\_height** (*A, B, C*)

Parameters are coordinates of points which are tops of triangle. The function calculates the ABC triangle height.

**Parameters**

- **A** – coordinates of the first point [A top]
- **B** – coordinates of the second point [B top]
- **C** – coordinates of the third point [C top]

**Returns** one value of ABC triangle height

**vectors\_angle** (*A, B*)

This function calculates the angle between two given vectors (starting from the [0,0,0] to the given coordinates).

**Parameters**

- **A** – coordinates of the first point which is the end of the vector
- **B** – coordinates of the second point which is the end of the vector

**Returns** the angle between vectors in question (in radians)

**vectors\_angle\_alt** (*A, B*)

This function calculates the angle between two given vectors (starting from the [0,0,0] to the given coordinates)

- alternative method.

**Parameters**

- **A** – coordinates of the first point which is the end of the vector
- **B** – coordinates of the second point which is the end of the vector

**Returns** the angle between vectors in question (in radians)

**vectors\_angle\_alt\_anorm** (*A, B, A\_norm*)

**This function calculates the angle between two given vectors (starting from the [0,0,0] to the given coordinates**

- alternative method with additional *A\_norm* holding norm of *A*.

**Parameters**

- **A** – coordinates of the first point which is the end of the vector
- **B** – coordinates of the second point which is the end of the vector
- **A\_norm** – additional parameter holding normalized of vector *A*

**Returns** the angle between vectors in question (in radians)

**vectors\_angle\_anorm** (*A, B, A\_norm*)

**This function calculates the angle between two given vectors (starting from the [0,0,0] to the given coordinates**  
using additional *A\_norm* holding norm of *A*.

**Parameters**

- **A** – coordinates of the first point which is the end of the vector
- **B** – coordinates of the second point which is the end of the vector
- **A\_norm** – additional parameter holding normalized of vector A

**Returns** the angle between vectors in question (in radians)

**class LinearizeOneWay**

Bases: `object`

**here** (*coords*)

This function simplifies the trace by removing the redundant, linear points :param coords: 3D coordinates of a trace as an array-like object :return: indices of coordinates which are a staring and ending points of linear fragments and other non-linear points of the trace

**class LinearizeHobbit**

Bases: `aquaduct.geom.traces.LinearizeOneWay`

**and\_back\_again** (*coords*)

**\_\_call\_\_** (...) <==> x(...)

**class LinearizeRecursive**

Bases: `object`

Base class for linearization methods classes.

It implements recursive algorithm.

**here** (*coords*, *depth=0*)

Core of recursive linearization algorithm.

It checks if the first, the last and the middle point are linear according to the criterion. The middle point is a selected point that is in the middle of length of the paths made by input coordinates.

If these points are linear their indices are returned. Otherwise, coordinates are split into two parts. First part spans points from the first point to the middle point (inclusive) and the second part spans points from the middle (inclusive) to the last point. Next, these two parts are submitted recursively to `here()`.

Results of these recursive calls are joined, redundant indices are removed and sorted result is returned.

**Parameters**

- **coords** (`numpy.ndarray`) – Input coordinates.
- **depth** (`int`) – Depth of recurrence.

**Returns** Indices of `coords` points that can be used instead of all points in visulatization.

**Return type** list of int

**\_\_call\_\_** (...) <==> x(...)

**class TriangleLinearize** (*threshold=0.01*)

Bases: `object`

**\_\_init\_\_** (*threshold=0.01*)

x.**\_\_init\_\_**(...) initializes x; see help(type(x)) for signature

**is\_linear** (*coords*, *\*\*kwargs*)

**class VectorLinearize** (*threshold=0.05236*)

Bases: `object`

Base class for linearization methods classes.

It implements vector linearization criterion.

`__init__` (*threshold=0.05236*)

`x.__init__`(...) initializes x; see `help(type(x))` for signature

**`is_linear_core`** (*coords, depth=None*)

Method checks if input coordinates are linear according to the threshold and depth.

It begins with calculation of the threshold. If *depth* is None it is set to 1. Current threshold is calculated with following simple equation:

$$threshold_{current} = threshold_{initial} * (2 - 0.9^{depth})$$

Next, in a loop over all points but the first and the last the angle is calculated between two vectors. The first one made by the point and the first point, and the second vector made by the last and the first point. If any of the calculated angles is bigger the the threshold methods returns False; otherwise method returns True.

#### Parameters

- **`coords`** (*numpy.ndarray*) – Coordinates for which linearization criterion is checked.
- **`depth`** (*int*) – Depth of recurrence.

**Returns** True if input coordinates are linear and False otherwise.

**Return type** `bool`

**`is_linear`** (*coords, depth=None, \*\*kwargs*)

For more detail see `is_linear_core()` which is used as the criterion of linearity in this method.

#### Parameters

- **`coords`** (*numpy.ndarray*) – Coordinates for which linearization criterion is checked.
- **`depth`** (*int*) – Depth of recurrence.

**Returns** True if input coordinates are linear and False otherwise. Criterion is checked for coordinates in normal and reverse order.

**Return type** `bool`

**`class LinearizeRecursiveVector`** (*threshold=0.05236*)

Bases: `aquaduct.geom.traces.LinearizeRecursive`, `aquaduct.geom.traces.VectorLinearize`

Class provides recursive linearization of coordinates with `LinearizeRecursive` algorithm and the criterion of linearity implemented by `VectorLinearize`. This is default method.

**`class LinearizeRecursiveTriangle`** (*threshold=0.01*)

Bases: `aquaduct.geom.traces.LinearizeRecursive`, `aquaduct.geom.traces.TriangleLinearize`

Class provides recursive linearization of coordinates with `LinearizeRecursive` algorithm and the criterion of linearity implemented by `TriangleLinearize`.

**`class LinearizeHobbitVector`** (*threshold=0.05236*)

Bases: `aquaduct.geom.traces.LinearizeHobbit`, `aquaduct.geom.traces.VectorLinearize`

Class provides recursive linearization of coordinates with `LinearizeHobbit` algorithm and the criterion of linearity implemented by `VectorLinearize`.

**`class LinearizeHobbitTriangle`** (*threshold=0.01*)

Bases: `aquaduct.geom.traces.LinearizeHobbit`, `aquaduct.geom.traces.TriangleLinearize`

Class provides recursive linearization of coordinates with `LinearizeHobbit` algorithm and the criterion of linearity implemented by `TriangleLinearize`.

**class LinearizeOneWayVector** (*threshold=0.05236*)

Bases: `aquaduct.geom.traces.LinearizeOneWay`, `aquaduct.geom.traces.VectorLinearize`

Class provides recursive linearization of coordinates with `LinearizeOneWay` algorithm and the criterion of linearity implemented by `VectorLinearize`.

**class LinearizeOneWayTriangle** (*threshold=0.01*)

Bases: `aquaduct.geom.traces.LinearizeOneWay`, `aquaduct.geom.traces.TriangleLinearize`

Class provides recursive linearization of coordinates with `LinearizeOneWay` algorithm and the criterion of linearity implemented by `TriangleLinearize`.

## Module contents

### 5.1.1.3 aquaduct.traj package

#### Submodules

#### aquaduct.traj.barber module

Module implements AutoBarber generation of spheres.

**class Sphere**

Bases: `aquaduct.traj.barber.Sphere`

Simple sphere class.

**is\_point\_within** (*point*)

**is\_sphere\_within** (*sphere*)

**is\_sphere\_cloud** (*sphere*)

**class WhereToCut** (*spaths=None*, *inlets=None*, *expected\_nr\_of\_spaths=None*, *selection=None*, *mincut=None*, *mincut\_level=False*, *maxcut=None*, *maxcut\_level=False*, *tovdw=False*, *forceempty=False*)

Bases: `aquaduct.traj.sandwich.ReaderAccess`

Class implements method for creating (optimal) set of AutoBarber spheres for a collection of spaths; access to trajectory is also required to read VdW radii.

**\_\_init\_\_** (*spaths=None*, *inlets=None*, *expected\_nr\_of\_spaths=None*, *selection=None*, *mincut=None*, *mincut\_level=False*, *maxcut=None*, *maxcut\_level=False*, *tovdw=False*, *forceempty=False*)

#### Parameters

- **spaths** (*list*) – List of `aquaduct.traj.paths.SinglePath` objects.
- **expected\_nr\_of\_spaths** (*int*) – Number of spaths passed. Required when length of spaths cannot be calculated, eg when it is a generator.
- **selection** (*str*) – Selection string of molecular object used for spheres generation.
- **mincut** (*float*) – Value of *mincut* parameter.
- **maxcut** (*float*) – Value of *maxcut* parameter.
- **mincut\_level** (*bool*) – Flag of *mincut\_level*.
- **maxcut\_level** (*bool*) – Flag of *maxcut\_level*.
- **tovdw** (*bool*) – Flag of to VdW radii correction parameter.



- **forceempty** (*bool*) – If set *True* spheres of radius 0 are returned if no other sphere can be generated.

```

check_minmaxcuts ()
add_spheres_from_spaths (spaths)
add_spheres_from_inlets (inlets)
get_current_nr ()
inlet2sphere (inlet)
spath2spheres (sp)
_cut_thyself (spheres_passed, progress=False)
cut_thyself ()
is_overlapping_with_cloud (sphere)
cloud_groups (progress=False)

```

### aquaduct.traj.dumps module

```

class TmpDumpWriterOfMDA
    Bases: object
    __init__ ()
        x.__init__(...) initializes x; see help(type(x)) for signature
    dump_frames (reader, frames, selection='protein')
    close ()
    __del__ ()

```

### aquaduct.traj.inlets module

```

class ProtoInletTypeCodes

    surface = 'surface'
    internal = 'internal'
    incoming = 'inin'
    outgoing = 'inout'

class InletTypeCodes
    Bases: aquaduct.traj.inlets.ProtoInletTypeCodes
    all_surface = [('surface', 'inin'), ('surface', 'inout')]
    all_internal = [('internal', 'inin'), ('internal', 'inout')]
    all_incoming = [('surface', 'inin'), ('internal', 'inin')]
    all_outgoing = [('surface', 'inout'), ('internal', 'inout')]
    surface_incoming = ('surface', 'inin')
    internal_incoming = ('internal', 'inin')
    internal_outgoing = ('internal', 'inout')
    surface_outgoing = ('surface', 'inout')
    itype = 'internal'

```

```
class InletClusterGenericType (inp, out)
    Bases: object
    __init__ (inp, out)
        x.__init__(...) initializes x; see help(type(x)) for signature
    input
    output
    static cluster2str (cl)
    __getitem__ (item)
    __len__ ()
    __str__ () <==> str(x)
    __repr__ () <==> repr(x)
    make_val (base)
    __cmp__ (other)
    __hash__ () <==> hash(x)

class InletClusterExtendedType (surfin, interin, interout, surfout)
    Bases: aquaduct.traj.inlets.InletClusterGenericType
    __init__ (surfin, interin, interout, surfout)
        x.__init__(...) initializes x; see help(type(x)) for signature
    generic

class Inlet (coords=None, type=None, reference=None, frame=None)
    Bases: object
    __init__ (coords=None, type=None, reference=None, frame=None)
        x.__init__(...) initializes x; see help(type(x)) for signature

class Inlets (spath, center_of_system=None, onlytype=[('surface', 'inin'), ('surface', 'inout')],
    passing=False, pbar=None)
    Bases: object
    __init__ (spath, center_of_system=None, onlytype=[('surface', 'inin'), ('surface', 'inout')],
    passing=False, pbar=None)
        x.__init__(...) initializes x; see help(type(x)) for signature
    add_leaf_wrapper (name=None, message=None, toleaf=None)
    resize_leaf_0 ()
    add_message_wrapper (message=None, toleaf=None)
    extend_inlets (spath, onlytype=None)
    add_cluster_annotations (clusters)
    add_outliers_annotations (new_clusters)
    add_spheres (spheres)
    get_inlets_references ()
    size
    coords
    types
    refs
    refs_names
```

```
call_clusterization_method (method, data, spheres=None)  
get_flat_tree (message=None)  
perform_clustering (method)  
perform_reclustering (method, skip_outliers=False, skip_size=None)  
recluster_cluster (method, cluster)  
recluster_outliers (method)  
small_clusters_to_outliers (maxsize)  
renumber_clusters ()  
sort_clusters ()  
clusters_list  
clusters_centers  
clusters_size  
clusters_std  
spath2ctypes (**kwargs)  
spath2ctype (sp)  
lim_to (what, towhat)  
lim2spaths (spaths)  
lim2rnames (rnames)  
lim2types (types)  
lim2clusters (clusters)  
limspaths2 (**kwargs)  
get_chull ()
```

### **aquaduct.traj.paths module**

```
union_full (a, b)  
union_smarter (a, b)  
union (a, b, smarter=True)  
glue (a, b)  
xor_full (*args, **kwargs)  
xor_smarter (*args, **kwargs)  
xor (a, b, smarter=True)  
left (a, b, smarter=True)  
right (a, b, smarter=True)  
class PathTypesCodes  
  
    path_in_code = 'i'  
    path_object_code = 'c'  
    path_out_code = 'o'  
    path_walk_code = 'w'
```

```
class GenericPathTypeCodes
```

```
    object_name = 'c'
```

```
    scope_name = 's'
```

```
    out_name = 'n'
```

```
class GenericPaths (id_of_res, name_of_res=None, single_res_selection=None, min_pf=None,
                    max_pf=None)
```

```
    Bases: object, aquaduct.traj.paths.GenericPathTypeCodes
```

```
    __init__ (id_of_res, name_of_res=None, single_res_selection=None, min_pf=None,
              max_pf=None)
```

```
        x.__init__(...) initializes x; see help(type(x)) for signature
```

```
    types
```

```
    frames
```

```
    coords
```

```
    max_frame
```

```
    min_frame
```

```
    add_012 (os_in_frames)
```

```
    add_object (frame)
```

```
    add_scope (frame)
```

```
    add_type (frame, ftype)
```

```
    _gpt ()
```

```
    _gpo ()
```

```
    _gpi ()
```

```
    get_paths_in ()
```

```
    get_paths_out ()
```

```
    find_paths (fullonly=False, smartr=True)
```

```
    find_paths_types (fullonly=False)
```

```
    get_single_path_types (spath)
```

```
    barber_with_spheres (spheres)
```

```
class SinglePathID (path_id=None, nr=None, name=None)
```

```
    Bases: object
```

```
    __init__ (path_id=None, nr=None, name=None)
```

```
        x.__init__(...) initializes x; see help(type(x)) for signature
```

```
    __str__ () <==> str(x)
```

```
    __eq__ (other)
```

```
        x.__eq__(y) <==> x==y
```

```
yield_single_paths (gps, fullonly=None, progress=None, passing=None)
```

```
yield_generic_paths (spaths, progress=None)
```

```
class MacroMolPath (path_id, paths, types, single_res_selection=None)
```

```
    Bases: object, aquaduct.traj.paths.PathTypesCodes, aquaduct.traj.inlets.InletTypeCodes
```

```
    empty_coords = array([], shape=(0, 3), dtype=float64)
```

`__init__(path_id, paths, types, single_res_selection=None)`  
x.`__init__(...)` initializes x; see `help(type(x))` for signature

`object_len`

`is_single()`

`is_passing()`

`is_frame_in(frame)`

`is_frame_object(frame)`

`is_frame_out(frame)`

`is_frame_walk(frame)`

`path_in`

`path_object`

`path_out`

`types_in`

`types_object`

`types_out`

`coords_first_in`

`paths_first_in`

`coords_last_out`

`paths_last_out`

`coords_filo`

`get_inlets()`

`coords_in`

`coords_object`

`coords_out`

`coords`

`coords_cont`

`_paths`

`paths`

`paths_cont`

`types`

`types_cont`

`gtypes`

`gtypes_cont`

`etypes`

`etypes_cont`

`size`

`sizes`

`begins`

`ends`

```
has_in
has_object
has_out
get_coords (**kwargs)
_make_smooth_coords (smooth)
get_coords_cont (smooth=None)
get_distance_cont (smooth=None, normalize=False)
get_distance_rev_cont (*args, **kwargs)
get_distance_both_cont (**kwargs)
get_velocity_cont (**kwargs)
get_acceleration_cont (**kwargs)
_MacroMolPath__object_len_calculate ()
class SinglePath (path_id, paths, types, single_res_selection=None)
    Bases: aquaduct.traj.paths.MacroMolPath
    is_single ()
    is_passing ()
class PassingPath (path_id, paths, types, single_res_selection=None)
    Bases: aquaduct.traj.paths.MacroMolPath
    __init__ (path_id, paths, types, single_res_selection=None)
        x.__init__(...) initializes x; see help(type(x)) for signature
    __has_out = None
        self.has_in = True self.has_out = True
    object_len
    has_in
    has_out
    is_single ()
    is_passing ()
    is_frame_walk (frame)
    types
    gtypes
    sizes
    _paths
    coords
    path
    paths
    coords_first_in
    paths_first_in
    coords_last_out
    paths_last_out
    get_coords (smooth=None)
```

```

    get_inlets ()
    _PassingPath__object_len_calculate ()
class MasterPath (sp)
    Bases: aquaduct.traj.paths.MacroMolPath
    __init__ (sp)
        x.__init__(...) initializes x; see help(type(x)) for signature
    add_width (width)

```

### aquaduct.traj.sandwich module

```

class Window (start, stop, step)
    Bases: object
    __init__ (start, stop, step)
        x.__init__(...) initializes x; see help(type(x)) for signature
    __repr__ () <==> repr(x)
    range ()
    get_real (frame)
    len ()

```

```

class MasterReader
    Bases: object
    open_reader_traj = {}
    topology = ''
    trajectory = []
    window = None
    sandwich_mode = None
    engine_name = 'mda'
    __call__ (topology, trajectory, window=None, sandwich=False)

```

#### Parameters

- **topology** (*str*) – Topology file name.
- **trajectory** (*list*) – List of trajectories. Each element is a fine name.
- **window** (*Window*) – Frames window to read.
- **sandwich** (*bool*) – Flag for setting sandwich mode.

```

__getstate__ ()
__setstate__ (state)
engine
correct_window ()
__repr__ () <==> repr(x)
sandwich (number=False)
baguette (number=False)
iterate (number=False)
get_single_reader (number)

```

```
get_reader_by_id(someid)
real_number_of_frames()
number_of_frames(one_layer=False)
number_of_layers()
```

```
class ReaderAccess
```

```
    Bases: object
```

```
    reader
```

```
VdW_radii = {'ac': 2.47, 'ag': 2.11, 'al': 1.84, 'am': 2.44, 'ar': 1.88, 'as': 1.8}
    Dictionary of VdW radii.
```

Data taken from L. M. Mentel, mendelev, 2014. Available at: <https://bitbucket.org/lukaszmentel/mendelev>. Package **mendelev** is not used because it depends on too many other libraries.

```
class ReaderTraj(topology, trajectory, number=None, window=None)
```

```
    Bases: aquaduct.traj.sandwich.ReaderAccess
```

```
    __init__(topology, trajectory, number=None, window=None)
```

#### Parameters

- **topology** (*str*) – Topology file name.
- **trajectory** (*list*) – Trajectory file name.
- **number** (*int*) – Number of trajectory file.
- **window** (*Window*) – Frames window to read.
- **reader** (*Reader*) – Parent reader object.

This is base class for MD data access engines.

```
__repr__() <==> repr(x)
```

```
iterate_over_frames()
```

```
set_frame(frame)
```

```
dump_frames(frames, selection=None, filename=None)
```

```
__del__()
```

```
open_trajectory()
```

```
close_trajectory()
```

```
set_real_frame(real_frame)
```

```
real_number_of_frames()
```

```
parse_selection(selection)
```

```
atom_vdw(atomid)
```

```
atom2residue(atomid)
```

```
atoms_positions(atomids)
```

```
residues_positions(resids)
```

```
residues_names(resids)
```

```
topology_resids(resids)
```

```
atoms_masses(atomids)
```

```
dump_frames_to_file(frames, selection, filename)
```

```
class ReaderTrajViaMDA(topology, trajectory, number=None, window=None)
```

```
    Bases: aquaduct.traj.sandwich.ReaderTraj
```



```
open_trajectory ()
close_trajectory ()
set_real_frame (real_frame)
parse_selection (selection)
atom2residue (atomid)
atoms_positions (atomids)
residues_positions (resids)
residues_names (resids)
topology_resids (resids)
real_number_of_frames ()
atoms_masses (atomids)
atom_vdw (atomid)
dump_frames_to_file (frames, selection, filename)

class Selection (selected)
    Bases: aquaduct.traj.sandwich.ReaderAccess
    __init__ (selected)
        x.__init__(...) initializes x; see help(type(x)) for signature
    layer (number)
    numbers ()
    ix (ix)
    len ()
    get_reader (number)
    add (other)
    uniquify ()
    ids ()
    coords ()
    center_of_mass ()

class AtomSelection (selected)
    Bases: aquaduct.traj.sandwich.Selection
    vdw ()
    residues ()
    coords ()
    center_of_mass ()
    contains_residues (other_residues, convex_hull=False, map_fun=None, known_true=None)
    containing_residues (other_residues, *args, **kwargs)
    chull ()

class ResidueSelection (selected)
    Bases: aquaduct.traj.sandwich.Selection
    coords ()
    names ()
```

```
single_residues ()
coords_range_core (*args, **kwargs)
coords_range (srange, number, rid)
class FramesRangeCollection
    Bases: object
    __init__ ()
        x.__init__(...) initializes x; see help(type(x)) for signature
    append (srange)
    get_ranges (srange)
smooth_coords_ranges (*args, **kwargs)
class SingleResidueSelection (resid)
    Bases: aquaduct.traj.sandwich.ReaderAccess
    __init__ (resid)
        x.__init__(...) initializes x; see help(type(x)) for signature
    topology_resid
    get_reader ()
    coords (frames)
    _coords (**kwargs)
    coords_smooth (sranges, smooth)
```

## Module contents

### 5.1.1.4 aquaduct.utils package

#### Submodules

#### aquaduct.utils.clui module

Module comprises conveniences functions and definitions for different operations related to command line user interface.

```
class roman_emulation
    Bases: object
    toRoman (nr)
emit_message_to_file_in_root_logger (mess)
message_special (mess)
message (mess, cont=False)
    Prints message to standard error. If FileHandler is present in the root_logger the same message is
    appended to the log file.
```

#### Parameters

- **mess** (*str*) – message to print
- **cont** (*bool*) – if set True no new line is printed

```
class fbm (info, cont=True)
    Bases: object
```

```

__init__(info, cont=True)
    x.__init__(...) initializes x; see help(type(x)) for signature
__enter__()
__exit__(typ, value, traceback)
__call__(...) <==> x(...)

```

**class tictoc** (*mess*)

Bases: `object`

```

__init__(mess)
    x.__init__(...) initializes x; see help(type(x)) for signature
__enter__()
__exit__(typ, value, traceback)

```

**gregorian\_year\_in\_days** = 365.2425

Length of Gregorian year in days. Average value. Source: <https://en.wikipedia.org/wiki/Year>

**smart\_time\_string** (*s, rl=0, t=1.1, maximal\_length=None, maximal\_units=5*)

Function transforms time in seconds to nicely formatted string of length defined by `maximal_length`. Depending on number of seconds time is represented with one or more of the following units:

Unit name	Unit abbreviation
seconds	s
minutes	m
hours	h
days	d
years	y

Maximal number of units used in time string can be set with `maximal_units`.

#### Parameters

- **s** (*int*) – Input time in seconds.
- **rl** (*int*) – Number of units already used for representing time.
- **t** (*float*) – Exces above standard number of current time units.
- **maximal\_length** (*int*) – Maximal length of the output string. Must be greater then 0.
- **maximal\_units** (*int*) – Maximal number of units used in the output string. Must be greater then 0 and lower then 6.

**Returns** string of nicely formatted time

**Return type** `str`

**gsep** (*sep='-', times=72, length=None*)

Generic separator.

#### Parameters

- **sep** (*str*) – Element(s) of separator.
- **times** (*int*) – Number of times `sep` is printed.
- **length** (*int*) – Optional maximal length of output.

**Returns** String separator.

**Return type** `str`

**tsep** (*line*)

**Parameters** `line (str)` – Input line.

**Returns** Returns default `gsep()` of length of `line`.

**underline (line)**

**Parameters** `line (str)` – Input line.

**Returns** String made by concatenation of `line`, `os.linesep`, and output of `tsep()` called with `line`.

**Return type** `str`

**thead (line)**

**Parameters** `line (str)` – Input line.

**Returns** String made by concatenation of output of `tsep()` called with `line`, `line`, `os.linesep`, and again output of `tsep()` called with `line`.

**Return type** `str`

**class SimpleProgressBar (maxval=None, mess=None)**

Bases: `object`

Simple progress bar displaying progress with percent indicator, progress bar and ETA. Progress is measured by iterations.

#### Variables

- `rotate (str)` – String comprising characters with frames of a rotating toy.
- `barlength (int)` – Length of progress bar.
- `maxval (int)` – maximal number of iterations
- `current (int)` – current number of iterations
- `overrun_notice (bool)` – if True, overrun above `maxval` iterations causes insert of newline
- `overrun (bool)` – flag of overrun
- `begin (int)` – time in seconds at the initialization of the `SimpleProgressBar` class.
- `tcurrent (int)` – time in seconds of current iteration

`rotate = '\\|/-'`

`barlength = 24`

`__init__ (maxval=None, mess=None)`

#### Parameters

- `maxval (int)` – Maximal number of iterations stored to `maxval`.
- `mess (str)` – Optional message displayed at progress bar initialization.

`bar ()`

`ETA ()`

Returns ETA calculated on the basis of current number of iterations `current` and current time `tcurrent`. If number of iterations is 0 returns `?`. Time is formatted with `smart_time_string()`.

**Returns** ETA as string.

**Return type** `str`

`percent ()`

Returns float number of percent progress calculated in the basis of current number of iterations `current`. Should return number between 0 and 100.

**Returns** percent progress number

**Return type** float

**show()**

Shows current progress.

If value returned by `percent()` is  $\leq 100$  then progress is printed as percent indicator leaded by ETA calculated by `ETA()`.

If value returned by `percent()` is  $> 100$  then progress is printed as number of iterations and total time.

Progress bar is written to standard error.

**heartbeat()**

**next()**

**update(step)**

Updates number of current iterations `current` by one if `step` is  $> 0$ . Otherwise number of current iterations is not updated. In both cases time of current iteration `tcurrent` is updated and `show()` is called.

**Parameters** `step(int)` – update step

**ttime()**

Calculates and returns total time string formatted with `smart_time_string()`.

**Returns** string of total time

**Return type** str

**finish()**

Finishes progress bar. First, `update()` is called with `step = 0`. Next message of total time is written to standard error.

**pbar**

alias of `aquaduct.utils.clui.SimpleProgressBar`

**get\_str\_timestamp()**

**class SimpleTree** (`name=None, message=None`)

Bases: `object`

**\_\_init\_\_(name=None, message=None)**

`x.__init__(...)` initializes `x`; see `help(type(x))` for signature

**\_\_repr\_\_()**  $\leq \Rightarrow$  `repr(x)`

**is\_leaf()**

**leafs\_names**

**get\_leaf(name)**

**add\_message** (`message=None, toleaf=None, replace=False`)

**add\_message\_to\_leaf** (`message=None, toleaf=None, replace=False`)

**add\_leaf** (`name=None, message=None, toleaf=None`)

**add\_leaf\_to\_leaf** (`name=None, message=None, toleaf=None`)

**print\_simple\_tree** (`st, prefix=None, multiple=False, concise=True`)

## aquaduct.utils.helpers module

Collection of helpers - functions and decorators.

**combine** (*seqin*)

This is an alien function. It is not extensively used.

Directly taken from [http://aspn.activestate.com/ASPN/Cookbook/Python/Recipe/302478/index\\_txt](http://aspn.activestate.com/ASPN/Cookbook/Python/Recipe/302478/index_txt)

Returns a list of all combinations of argument sequences. For example, following call:

```
combine(( (1,2), (3,4) ))
```

gives following list of combinations:

```
[[1, 3], [1, 4], [2, 3], [2, 4]]
```

**Parameters** *seqin* (*tuple*) – Tuple of sequences to combine.

**Returns** All possible combinations of all input sequences.

**Return type** list of lists

**are\_rows\_uniq** (*some\_array*)**robust\_and** (*a*, *b*)**robust\_or** (*a*, *b*)**is\_number** (*s*)**lind** (*l*, *ind*)

Indexes lists using lists of integers as identifiers. For example:

```
lind(['a', 'b', 'c', 'd', 'e'], [1,4,2])
```

returns:

```
['b', 'e', 'c']
```

**Parameters**

- **l** (*list*) – List to be indexed.
- **ind** (*list*) – Integer indexes.

**Returns** Reindexed list.

**Return type** list

**class Auto**

Auto type definition. The class is used as an alternative value for options (if particular option supports it). If options (or variables/parameters etc.) have value of *Auto* it means that an automatic process for parametrization should be performed.

For example, if the input parameter is set to *Auto* it is supposed that its value is calculated on the basis of input data or other parameters.

**\_\_repr\_\_** ()

**Returns** String *Auto*.

**Return type** *str*

**\_\_str\_\_** ()

Calls *\_\_repr\_\_* ().

**create\_tmpfile** (*ext=None*)

Creates temporary file. File is created, closed and its file name is returned.

---

**Note:** It is responsibility of the caller to delete the file.

---

**Parameters** `ext` (*str*) – Optional extension of the file.

**Returns** File name of created temporary file.

**Return type** *str*

**range2int** (*r*, *uniq=True*)

Transforms a string range in to a list of integers (with added missing elements from given ranges).

For example, a following string:

```
'0:2 4:5 7 9'
```

is transformed into:

```
[0, 1, 2, 4, 5, 7, 9]
```

**Parameters**

- `r` (*str*) – String of input range.
- `uniq` (*bool*) – Optional parameter, if set to *True* only unique and sorted integers are returned.

**Returns** List of integers.

**Return type** list of int

**int2range** (*l*)

Transforms a list of integers in to a string of ranges.

For example, a following list:

```
[0, 1, 2, 4, 5, 7, 9]
```

is transformed into:

```
0:2 4:5 7 9
```

**Parameters** `l` (*list*) – input list of int

**Returns** String of ranges.

**Return type** *str*

**is\_iterable** (*l*)

Checks if provided object is iterable. Returns True if it is iterable, otherwise returns False.

**Parameters** `l` (*list*) – input object

**Returns** True if submitted object is iterable otherwise returns False.

**Return type** *bool*

**sortify** (*gen*)

Decorator to convert functions' outputs into a sorted list. If the output is iterable it is converted in to a list of appropriate length. If the output is not iterable it is converted in to a list of length 1.

Written on the basis of `listify()`.

**Returns** Output of decorated function converted to a sorted list.

**Return type** list

**uniqify** (*gen*)

Decorator to convert functions' outputs into a sorted list of unique objects. If the output is iterable it is converted in to a list of appropriate length. If the output is not iterable it is converted in to a list of length 1.

Written on the basis of `listify()`.

**Returns** Output of decorated function converted to a sorted list of unique objects.

**Return type** list

**noaction** (*gen*)**listify** (*gen*)

Decorator to convert functions' outputs into a list. If the output is iterable it is converted in to a list of appropriate length. If the output is not iterable it is converted in to a list of length 1.

This function was copied from:

<http://argandgahandapandpa.wordpress.com/2009/03/29/python-generator-to-list-decorator/>

and further improved by `tljm@wp.pl`.

**Returns** Output of decorated function converted to a list.

**Return type** list

**tupleify** (*gen*)

Decorator to convert functions' outputs into a tuple. If the output is iterable it is converted in to a tuple of appropriate length. If the output is not iterable it is converted in to a tuple of length 1.

Written on the basis of `listify()`.

**Returns** Output of decorated function converted to a tuple.

**Return type** tuple

**class arrayify** (*shape=None*)

Bases: `object`

`__init__` (*shape=None*)

`x.__init__(...)` initializes x; see `help(type(x))` for signature

`__call__` (*gen*)

Decorator to convert functions' outputs into a 2D numpy array. If the output is iterable it is converted in to a 2D numpy array of appropriate shape. If the output is not iterable it is converted in to a 2D numpy array of shape 1x1.

Written on the basis of `listify()`.

**Returns** Output of decorated function converted to a 2D numpy array.

**Return type** `numpy.ndarray`

**arrayify1** (*gen*)

Decorator to convert functions' outputs into a 1D numpy array. If the output is iterable it is converted in to a 2D numpy array of appropriate shape. If the output is not iterable it is converted in to a 2D numpy array of shape 1x1.

Written on the basis of `listify()`.

**Returns** Output of decorated function converted to a 1D numpy array.

**Return type** `numpy.ndarray`

**list\_blocks\_to\_slices** (*l*)

Slices list in to block according to its elements identity. Resulting slices correspond to blocks of identical elements.

**Parameters** **1** (*list*) – List of any objects.

**Returns** Generator of slices.



**Return type** generator

**split\_list** (*l, s*)

**what2what** (*what, towhat*)

This function search if elements of the one list (:attr: 'what') are present in the other list (:attr: 'towhat') and returns indices of elements form :attr:'what' list as a tuple. If elements from the first list are not present in the second list the tuple is empty. :param list what: Input list for which indices of elements present in towhat are returned. :param list towhat: List of elements which input list is indexed to. :return: Indices of what list that are present in towhat list. :rtype: tuple

**make\_iterable** (*something*)

If input object is not iterable returns it as one element list. Otherwise returns the object.

**Parameters** *something* (*object*) – Input object.

**Returns** Iterable object.

**Return type** iterable or list

**iterate\_or\_die** (*something, times=None*)

**stretch\_zip** (*\*args*)

**compress\_zip** (*\*args*)

**zip\_zip** (*\*args, \*\*kwargs*)

**xzip\_xzip** (*\*args, \*\*kwargs*)

**concatenate** (*\*args*)

Concatenates input iterable arguments in to one generator.

**class Bunch** (*\*\*kws*)

Bases: *object*

<http://code.activestate.com/recipes/52308> foo=Bunch(a=1,b=2)

**\_\_init\_\_** (*\*\*kws*)

x.\_\_init\_\_(...) initializes x; see help(type(x)) for signature

**class SmartRangeFunction** (*element, times*)

Bases: *object*

**\_\_init\_\_** (*element, times*)

x.\_\_init\_\_(...) initializes x; see help(type(x)) for signature

**\_\_str\_\_** () <==> str(x)

**\_\_repr\_\_** () <==> repr(x)

**\_\_len\_\_** ()

**get** ()

**rev** ()

**isin** (*element*)

**first\_element** ()

**last\_element** ()

**overlaps** (*srange*)

**overlaps\_mutual** (*srange*)

**contains** (*srange*)

**class SmartRangeEqual** (*element, times*)

Bases: *aquaduct.utils.helpers.SmartRangeFunction*

**type** = 'e'

```
    get ()
    rev ()
    isin (element)
    last_element ()

class SmartRangeIncrement (element, times)
    Bases: aquaduct.utils.helpers.SmartRangeFunction
    type = 'i'
    get ()
    rev ()
    isin (element)
    last_element ()

class SmartRangeDecrement (element, times)
    Bases: aquaduct.utils.helpers.SmartRangeFunction
    type = 'd'
    get ()
    rev ()
    isin (element)
    last_element ()

class SmartRange (iterable=None)
    Bases: object
    __init__ (iterable=None)
        x.__init__(...) initializes x; see help(type(x)) for signature
    __str__ () <==> str(x)
    __repr__ () <==> repr(x)
    first_element ()
    last_element ()
    last_times ()
    raw
    append (element)
    get ()
    rev ()
    __len__ ()
    __iter__ ()
    min ()
    max ()
    isin (element)
```

## **aquaduct.utils.maths module**

**class** **NumpyDefaultsStorageTypes**

Bases: `object`

**float\_default**

alias of `numpy.float64`

**int\_default**

alias of `numpy.int64`

**make\_default\_array** (*array\_like*)

## **aquaduct.utils.multip module**

**class** **CpuThreadsCount**

Bases: `object`

**cpu\_count** = 2

**threads\_count** = None

## **Module contents**

### **5.1.1.5 aquaduct.visual package**

#### **Submodules**

#### **aquaduct.visual.cmaps module**

#### **aquaduct.visual.helpers module**

**euclidean** (*A, B*)

**cityblock** (*A, B*)

**cc\_safe** (*c*)

**cc** (*c*)

**color\_codes** (*code, custom\_codes=None*)

**get\_cmap** (*size*)

**class** **ColorMapDistMap**

Bases: `object`

**grey** = (0.5, 0.5, 0.5, 1)

**\_\_init\_\_** ()

`x.__init__(...)` initializes x; see `help(type(x))` for signature

**distance** (*E1, E2*)

**static color\_distance** (*e1, e2*)

**\_\_call\_\_** (...)  $\Leftrightarrow x(...)$

**\_ColorMapDistMap\_\_do\_cadex** ()

**f\_like** (*n*)

**aquaduct.visual.pymol\_cgo module****aquaduct.visual.pymol\_connector module**

```
class BasicPymolCGO
    Bases: object

    cgo_entity_begin = []
    cgo_entity_end = []

    __init__()
        x.__init__(...) initializes x; see help(type(x)) for signature

    clean (empty=False)

    new()

    get()

    static make_color_triple (color_definition)

class BasicPymolCGOLines
    Bases: aquaduct.visual.pymol_connector.BasicPymolCGO

    cgo_entity_begin = [2.0, 1.0]
    cgo_entity_end = [3.0]
    add (coords=None, color=None)

class BasicPymolCGOSpheres
    Bases: aquaduct.visual.pymol_connector.BasicPymolCGO

    cgo_entity_begin = []
    cgo_entity_end = []
    add (coords=None, radius=None, color=None)

class BasicPymolCGOPointers
    Bases: aquaduct.visual.pymol_connector.BasicPymolCGO

    cgo_entity_begin = []
    cgo_entity_end = []
    add_cone (coords1=None, coords2=None, radius1=None, radius2=None, color1=None, color2=None)
    add_pointer (point=None, direction=None, length=None, color=None, reverse=False)

class SimpleTarWriteHelper
    Bases: object

    __init__()
        x.__init__(...) initializes x; see help(type(x)) for signature

    open (filename)

    save_object2tar (obj, name)
    save_file2tar (filename, name)

    __del__()

class ConnectToPymol
    Bases: object

    cgo_line_width = 2.0
    ct_pymol = 'pymol'
```

```

ct_file = 'file'

__init__()
    x.__init__(...) initializes x; see help(type(x)) for signature

decode_color (**kwargs)

init_pymol()

init_script(filename)

add_cgo_object(name, cgo_object, state=None)

del_cgo_object(name, state=None)

load_pdb(name, filename, state=None)

orient_on(name)

__del__()

class SinglePathPlotter(pymol_connector, linearize=None)
    Bases: object

    __init__(pymol_connector, linearize=None)
        x.__init__(...) initializes x; see help(type(x)) for signature

    add_single_path_continuous_trace(spath, smooth=None, plot_in=True,
                                     plot_object=True, plot_out=True, plot_walk=True,
                                     **kwargs)

    paths_trace(spaths, smooth=None, name='paths', state=None, **kwargs)

    paths_inlets(spaths, smooth=None, color=None, plot_in=True, plot_out=True, name='in-out-
                 let', state=None, **kwargs)

    scatter(coords, radius=0.4, color='r', name='scatter', state=None)

    convexhull(chull, color='m', name='convexhull', state=None)

```

### aquaduct.visual.quickplot module

```

yield_spath_len_and_smooth_diff_in_types_slices(sp, smooth=None,
                                                  smooth_len=None,
                                                  smooth_diff=None,
                                                  types='etypes')

plot_colorful_lines(x, y, c, **kwargs)

spaths_spectra(spaths, **kwargs)

plot_spath_spectrum(sp, **kwargs)

spath_spectrum(sp, **kwargs)

showit(gen)

get_ax3d(fig, sub=111)

class SimpleTracePlotter
    Bases: object

    plot_line(coords, color, **kwargs)

    single_trace(coords, color='r', **kwargs)

    path_trace(path, color=('r', 'g', 'b'), plot_in=True, plot_object=True, plot_out=True,
                **kwargs)

class SimpleProteinPlotter
    Bases: aquaduct.visual.quickplot.SimpleTracePlotter

```

```
protein_trace (protein, smooth=None, color=('c', 'm', 'y'), **kwargs)

class SimplePathPlotter
    Bases: aquaduct.visual.quickplot.SimpleTracePlotter

    single_path_traces (spaths, smooth=None, color=('r', 'g', 'b'), **kwargs)

class MPLTracePlotter
    Bases: aquaduct.visual.quickplot.SimplePathPlotter, aquaduct.visual.quickplot.SimpleProteinPlotter

    init_ax (**kwargs)

    plot_line (**kwargs)

    scatter (**kwargs)
```

## Module contents

### 5.1.2 Module contents

Aqua-Duct - a collection of tools to trace residues in MD simulation.

**version()**

Returns *aquaduct* version number.

**Returns** 3 element tuple of int numbers

**Return type** tuple

**version\_nice()**

Returns *aquaduct* version number as nicely formatted string.

**Returns** string composed on the basis of the number returned by *version()*.

**Return type** str

**greetings()**

Returns fancy greetings of *aquaduct*. It has a form of ASCII-like graphic. Currently it returns following string:

```
-----
~ ~ ~ A Q U A - D U C T ~ ~ ~
#####
####          #####          #####          ####
##           #####          #####          ##
#            ##           ##           #
#            ##           ##           #
#            ##           ##           #
#            ##           ##           #
-----
```

**Returns** *aquaduct* fancy greetings.

**Return type** str

## AQUA-DUCT CHANGELOG

- **0.5.14 (26.09.2018)**
  - Fixed problems with calculation convex hull shapes for empty selections.
- **0.5.9 (12.03.2018)**
  - Rewritten module for MD data access. Sandwich mode added.
  - Coordinates can be stored in cache directory, in memory or generated on demand.
  - Support for long trajectories.
  - Passing through paths are supported.
  - Improvements in visualization script.
  - Coordinates of residues are calculated as center of geometry.
  - Recommended MDAnalysis is set to  $\geq 0.16$  and  $< 0.17$ . Version 0.17 is supported but not recommended.
  - Bug fixes and code cleanup.
- **0.4.0 - 0.4.14 (20.11.2017) unofficial**
  - Uses newest MDAnalysis (0.16.2).
  - Steady improvement of documentation (including API).
  - Names of traced molecules are returned in the result file and tables are split appropriately.
  - Tables in the result file are split in regard to Object and Passing paths.
  - Passing through paths are being introduced, WIP.
  - Additional tables in the result file.
  - CRD is enabled as topology/trajectory format.
  - Traced residues are identified by resindices instead of resids; this allows to use weak topologies such as PDB.
  - Removed roman dependency.
  - In addition to histograms approximate (ConvexHull approximation) areas and volumes of the scope and object can be calculated.
  - Bug fixes and reliability fixes.
- **0.3.7 (18.07.2017)**
  - Enable XTC trajectory format.
  - Reliability fix in progress bar display.
- **0.3.6 (28.06.2017)**
  - AQ can be run for given part of trajectory.

- Fixed bug in passing options to Barber clusterization method.
  - Recursive threshold can be defined as range; no disjoint ranges are supported.
- **0.3.5 (18.04.2017)**
  - As for now, the only supported version of MDAnalysis is 0.15.
- **0.3.4 (14.04.2017)**
  - Fixed bug in progress bar updating method causing critical error in some specific circumstances.
- **0.3.3 (20.03.2017)**
  - AutoBarber default values of maxcut\_level and mincut\_level changed to True.
  - Improved template configuration file.
  - Number of small improvements in documentation.
- **0.3.2 (24.02.2017)**
  - Major improvement: new auto\_barber based clustering method.
  - Clusterization history displayed as simple ascii tree.
  - AutoBarber min and max cut level options added.
  - Barber moved to separate module.
  - Fixed bug in visualization script; if no molecule is kept do not set style and color.
- **0.3.1 (04.02.2017)**
  - AutoBarber tovdw option.
  - AutoBarber minimal and maximal cut options.
  - Fixed bug in AutoBarber: some areas were sometimes not cut.
  - Documentation improvements.
  - Valve driver simplified. Most of the functionality moved to separate module.
  - Option for single precision storage.
  - Added Savitzky-Golay smoothing; AQ requires SciPy >= 0.14 now.
  - Improved sorting of CTypes.
  - Raw and Separate paths uses SmartRanges. This allowed for excellent performance improvement of Separate paths calculation.
  - Default display of molecule changed to silver cartoon.
  - Object shape displayed in orange.
  - Fixed several small bugs.
- **0.2.26 (21.01.2017)**
  - Stage execution time debug messages.
  - Total execution time debug message.
- **0.2.25 (18.01.2017)**
  - initial public release

PDF version of this documentation is also [available](#).

Go back to [current version](#) documentation.



## PYTHON MODULE INDEX

### a

- aquaduct, 74
- aquaduct.apps, 39
- aquaduct.apps.data, 37
- aquaduct.apps.valvecore, 38
- aquaduct.geom, 52
- aquaduct.geom.cluster, 39
- aquaduct.geom.convexhull, 40
- aquaduct.geom.master, 41
- aquaduct.geom.pca, 43
- aquaduct.geom.smooth, 44
- aquaduct.geom.traces, 48
- aquaduct.traj, 62
- aquaduct.traj.barber, 52
- aquaduct.traj.dumps, 53
- aquaduct.traj.inlets, 53
- aquaduct.traj.paths, 55
- aquaduct.traj.sandwich, 59
- aquaduct.utils, 71
- aquaduct.utils.clui, 62
- aquaduct.utils.helpers, 65
- aquaduct.utils.maths, 71
- aquaduct.utils.multip, 71
- aquaduct.visual, 74
- aquaduct.visual.cmaps, 71
- aquaduct.visual.helpers, 71
- aquaduct.visual.pymol\_cgo, 72
- aquaduct.visual.pymol\_connector, 72
- aquaduct.visual.quickplot, 73



## Symbols

- `_ColorMapDistMap__do_cadex()` (`ColorMapDistMap` method), 71
- `_MacroMolPath__object_len_calculate()` (`MacroMolPath` method), 58
- `_PassingPath__object_len_calculate()` (`PassingPath` method), 59
- `_ValveConfig__make_options_nt()` (`ValveConfig` method), 39
- `__call__()` (`CTypeSpathsCollectionWorker` method), 41
- `__call__()` (`Center` method), 43
- `__call__()` (`ColorMapDistMap` method), 71
- `__call__()` (`LinearizeHobbit` method), 50
- `__call__()` (`LinearizeRecursive` method), 50
- `__call__()` (`MasterReader` method), 59
- `__call__()` (`Normalize` method), 44
- `__call__()` (`PCA` method), 44
- `__call__()` (`PerformClustering` method), 40
- `__call__()` (`Smooth` method), 45
- `__call__()` (`Standartize` method), 44
- `__call__()` (`arrayify` method), 68
- `__call__()` (`fbm` method), 63
- `__cmp__()` (`InletClusterGenericType` method), 54
- `__del__()` (`ConnectToPymol` method), 73
- `__del__()` (`ReaderTraj` method), 60
- `__del__()` (`SimpleTarWriteHelper` method), 72
- `__del__()` (`TmpDumpWriterOfMDA` method), 53
- `__del__()` (`ValveDataAccessRoots` method), 38
- `__enter__()` (`fbm` method), 63
- `__enter__()` (`tictoc` method), 63
- `__eq__()` (`SinglePathID` method), 56
- `__exit__()` (`fbm` method), 63
- `__exit__()` (`tictoc` method), 63
- `__getitem__()` (`InletClusterGenericType` method), 54
- `__getstate__()` (`MasterReader` method), 59
- `__has_out` (`PassingPath` attribute), 58
- `__hash__()` (`InletClusterGenericType` method), 54
- `__init__()` (`ActiveWindowOverMaxStepSmooth` method), 47
- `__init__()` (`ActiveWindowSmooth` method), 46
- `__init__()` (`BarberClusterResult` method), 40
- `__init__()` (`BasicPymolCGO` method), 72
- `__init__()` (`Bunch` method), 69
- `__init__()` (`CTypeSpathsCollection` method), 42
- `__init__()` (`CTypeSpathsCollectionWorker` method), 41
- `__init__()` (`Center` method), 43
- `__init__()` (`ColorMapDistMap` method), 71
- `__init__()` (`ConnectToPymol` method), 73
- `__init__()` (`DistanceWindowOverMaxStepSmooth` method), 47
- `__init__()` (`DistanceWindowSmooth` method), 46
- `__init__()` (`FakeSingleResidueSelection` method), 43
- `__init__()` (`FloatWindow` method), 45
- `__init__()` (`FramesRangeCollection` method), 62
- `__init__()` (`GeneralWindow` method), 45
- `__init__()` (`GenericPaths` method), 56
- `__init__()` (`Inlet` method), 54
- `__init__()` (`InletClusterExtendedType` method), 54
- `__init__()` (`InletClusterGenericType` method), 54
- `__init__()` (`Inlets` method), 54
- `__init__()` (`IntWindow` method), 45
- `__init__()` (`LoadDumpWrapper` method), 37
- `__init__()` (`MacroMolPath` method), 56
- `__init__()` (`MasterPath` method), 59
- `__init__()` (`MaxStepSmooth` method), 47
- `__init__()` (`Normalize` method), 44
- `__init__()` (`PCA` method), 44
- `__init__()` (`PassingPath` method), 58
- `__init__()` (`PerformClustering` method), 40
- `__init__()` (`ReaderTraj` method), 60
- `__init__()` (`SavgolSmooth` method), 47
- `__init__()` (`Selection` method), 61
- `__init__()` (`SimpleProgressBar` method), 64
- `__init__()` (`SimpleTarWriteHelper` method), 72
- `__init__()` (`SimpleTree` method), 65
- `__init__()` (`SinglePathID` method), 56
- `__init__()` (`SinglePathPlotter` method), 73
- `__init__()` (`SingleResidueSelection` method), 62
- `__init__()` (`SmartRange` method), 70
- `__init__()` (`SmartRangeFunction` method), 69
- `__init__()` (`Smooth` method), 44
- `__init__()` (`Standartize` method), 44
- `__init__()` (`TmpDumpWriterOfMDA` method), 53
- `__init__()` (`TriangleLinearize` method), 50
- `__init__()` (`ValveConfig` method), 38
- `__init__()` (`ValveDataAccess_nc` method), 38
- `__init__()` (`VectorLinearize` method), 50
- `__init__()` (`WhereToCut` method), 52
- `__init__()` (`Window` method), 59
- `__init__()` (`WindowOverMaxStepSmooth` method), 47
- `__init__()` (`WindowSmooth` method), 46
- `__init__()` (`arrayify` method), 68

\_\_init\_\_() (fbm method), 62  
 \_\_init\_\_() (tictoc method), 63  
 \_\_iter\_\_() (SmartRange method), 70  
 \_\_len\_\_() (InletClusterGenericType method), 54  
 \_\_len\_\_() (SmartRange method), 70  
 \_\_len\_\_() (SmartRangeFunction method), 69  
 \_\_repr\_\_() (Auto method), 66  
 \_\_repr\_\_() (InletClusterGenericType method), 54  
 \_\_repr\_\_() (MasterReader method), 59  
 \_\_repr\_\_() (ReaderTraj method), 60  
 \_\_repr\_\_() (SimpleTree method), 65  
 \_\_repr\_\_() (SmartRange method), 70  
 \_\_repr\_\_() (SmartRangeFunction method), 69  
 \_\_repr\_\_() (Window method), 59  
 \_\_setstate\_\_() (MasterReader method), 59  
 \_\_str\_\_() (Auto method), 66  
 \_\_str\_\_() (InletClusterGenericType method), 54  
 \_\_str\_\_() (PerformClustering method), 40  
 \_\_str\_\_() (SinglePathID method), 56  
 \_\_str\_\_() (SmartRange method), 70  
 \_\_str\_\_() (SmartRangeFunction method), 69  
 \_coords() (SingleResidueSelection method), 62  
 \_cut\_thyself() (WhereToCut method), 53  
 \_edges() (in module aquaduct.geom.convexhull), 40  
 \_facets() (in module aquaduct.geom.convexhull), 40  
 \_get\_noclusters() (PerformClustering method), 40  
 \_gpi() (GenericPaths method), 56  
 \_gpo() (GenericPaths method), 56  
 \_gpt() (GenericPaths method), 56  
 \_make\_smooth\_coords() (MacroMolPath method), 58  
 \_paths (MacroMolPath attribute), 57  
 \_paths (PassingPath attribute), 58  
 \_point\_within\_convexhull() (in module aquaduct.geom.convexhull), 40  
 \_vertices\_ids() (in module aquaduct.geom.convexhull), 40  
 \_vertices\_points() (in module aquaduct.geom.convexhull), 40

## A

ActiveWindowOverMaxStepSmooth (class in aquaduct.geom.smooth), 47  
 ActiveWindowSmooth (class in aquaduct.geom.smooth), 46  
 add() (BasicPymolCGOLines method), 72  
 add() (BasicPymolCGOSpheres method), 72  
 add() (Selection method), 61  
 add\_012() (GenericPaths method), 56  
 add\_cgo\_object() (ConnectToPymol method), 73  
 add\_cluster\_annotations() (Inlets method), 54  
 add\_cone() (BasicPymolCGOPointers method), 72  
 add\_leaf() (SimpleTree method), 65  
 add\_leaf\_to\_leaf() (SimpleTree method), 65  
 add\_leaf\_wrapper() (Inlets method), 54  
 add\_message() (SimpleTree method), 65  
 add\_message\_to\_leaf() (SimpleTree method), 65  
 add\_message\_wrapper() (Inlets method), 54  
 add\_object() (GenericPaths method), 56

add\_outliers\_annotations() (Inlets method), 54  
 add\_pointer() (BasicPymolCGOPointers method), 72  
 add\_scope() (GenericPaths method), 56  
 add\_single\_path\_continuous\_trace() (SinglePathPlotter method), 73  
 add\_spheres() (Inlets method), 54  
 add\_spheres\_from\_inlets() (WhereToCut method), 53  
 add\_spheres\_from\_spaths() (WhereToCut method), 53  
 add\_type() (GenericPaths method), 56  
 add\_width() (MasterPath method), 59  
 all\_incoming (InletTypeCodes attribute), 53  
 all\_internal (InletTypeCodes attribute), 53  
 all\_outgoing (InletTypeCodes attribute), 53  
 all\_surface (InletTypeCodes attribute), 53  
 and\_back\_again() (LinearizeHobbit method), 50  
 append() (FramesRangeCollection method), 62  
 append() (SmartRange method), 70  
 aquaduct (module), 74  
 aquaduct.apps (module), 39  
 aquaduct.apps.data (module), 37  
 aquaduct.apps.valvecore (module), 38  
 aquaduct.geom (module), 52  
 aquaduct.geom.cluster (module), 39  
 aquaduct.geom.convexhull (module), 40  
 aquaduct.geom.master (module), 41  
 aquaduct.geom.pca (module), 43  
 aquaduct.geom.smooth (module), 44  
 aquaduct.geom.traces (module), 48  
 aquaduct.traj (module), 62  
 aquaduct.traj.barber (module), 52  
 aquaduct.traj.dumps (module), 53  
 aquaduct.traj.inlets (module), 53  
 aquaduct.traj.paths (module), 55  
 aquaduct.traj.sandwich (module), 59  
 aquaduct.utils (module), 71  
 aquaduct.utils.clui (module), 62  
 aquaduct.utils.helpers (module), 65  
 aquaduct.utils.maths (module), 71  
 aquaduct.utils.multip (module), 71  
 aquaduct.visual (module), 74  
 aquaduct.visual.cmaps (module), 71  
 aquaduct.visual.helpers (module), 71  
 aquaduct.visual.pymol\_cgo (module), 72  
 aquaduct.visual.pymol\_connector (module), 72  
 aquaduct.visual.quickplot (module), 73  
 aquaduct\_version\_nice() (in module aquaduct.apps.valvecore), 39  
 are\_rows\_uniq() (in module aquaduct.utils.helpers), 66  
 arrayify (class in aquaduct.utils.helpers), 68  
 arrayify1() (in module aquaduct.utils.helpers), 68  
 arrays2dict() (IdsOverIds static method), 38  
 atom2residue() (ReaderTraj method), 60  
 atom2residue() (ReaderTrajViaMDA method), 61  
 atom\_vdw() (ReaderTraj method), 60  
 atom\_vdw() (ReaderTrajViaMDA method), 61  
 atoms\_masses() (ReaderTraj method), 60  
 atoms\_masses() (ReaderTrajViaMDA method), 61  
 atoms\_positions() (ReaderTraj method), 60

atoms\_positions() (ReaderTrajViaMDA method), 61  
 AtomSelection (class in aquaduct.traj.sandwich), 61  
 Auto (class in aquaduct.utils.helpers), 66

## B

baguette() (MasterReader method), 59  
 bar() (SimpleProgressBar method), 64  
 barber\_with\_spheres() (GenericPaths method), 56  
 BarberCluster (class in aquaduct.geom.cluster), 40  
 BarberClusterResult (class in aquaduct.geom.cluster), 40  
 barlenght (SimpleProgressBar attribute), 64  
 BasicPymolCGO (class in aquaduct.visual.pymol\_connector), 72  
 BasicPymolCGOLines (class in aquaduct.visual.pymol\_connector), 72  
 BasicPymolCGOPointers (class in aquaduct.visual.pymol\_connector), 72  
 BasicPymolCGOSpheres (class in aquaduct.visual.pymol\_connector), 72  
 beat() (CTypeSpathsCollection method), 42  
 begins (MacroMolPath attribute), 57  
 Bunch (class in aquaduct.utils.helpers), 69

## C

cache (CoordsRangeIndexCache attribute), 37  
 cachedir (GlobalConfigStore attribute), 37  
 cachemem (GlobalConfigStore attribute), 37  
 call\_clusterization\_method() (Inlets method), 54  
 cc() (in module aquaduct.visual.helpers), 71  
 cc\_safe() (in module aquaduct.visual.helpers), 71  
 Center (class in aquaduct.geom.pca), 43  
 center\_of\_mass() (AtomSelection method), 61  
 center\_of\_mass() (Selection method), 61  
 centers() (PerformClustering method), 40  
 cgo\_entity\_begin (BasicPymolCGO attribute), 72  
 cgo\_entity\_begin (BasicPymolCGOLines attribute), 72  
 cgo\_entity\_begin (BasicPymolCGOPointers attribute), 72  
 cgo\_entity\_begin (BasicPymolCGOSpheres attribute), 72  
 cgo\_entity\_end (BasicPymolCGO attribute), 72  
 cgo\_entity\_end (BasicPymolCGOLines attribute), 72  
 cgo\_entity\_end (BasicPymolCGOPointers attribute), 72  
 cgo\_entity\_end (BasicPymolCGOSpheres attribute), 72  
 cgo\_line\_width (ConnectToPymol attribute), 72  
 check\_bounds\_at\_max\_window\_at\_pos() (GeneralWindow method), 45  
 check\_minmaxcuts() (WhereToCut method), 53  
 check\_version\_compliance() (in module aquaduct.apps.data), 37  
 check\_versions() (in module aquaduct.apps.data), 37  
 chull() (AtomSelection method), 61  
 cityblock() (in module aquaduct.visual.helpers), 71  
 clean() (BasicPymolCGO method), 72  
 close() (TmpDumpWriterOfMDA method), 53  
 close() (ValveDataAccess\_pickle method), 38

close\_all() (ValveDataAccessRoots method), 38  
 close\_trajectory() (ReaderTraj method), 60  
 close\_trajectory() (ReaderTrajViaMDA method), 61  
 cloud\_groups() (WhereToCut method), 53  
 cluster2str() (InletClusterGenericType static method), 54  
 cluster\_name() (ValveConfig static method), 38  
 clusters\_centers (Inlets attribute), 55  
 clusters\_list (Inlets attribute), 55  
 clusters\_size (Inlets attribute), 55  
 clusters\_std (Inlets attribute), 55  
 color\_codes() (in module aquaduct.visual.helpers), 71  
 color\_distance() (ColorMapDistMap static method), 71  
 ColorMapDistMap (class in aquaduct.visual.helpers), 71  
 combine() (in module aquaduct.utils.helpers), 65  
 common\_config\_names() (ValveConfig static method), 38  
 common\_traj\_data\_config\_names() (ValveConfig static method), 38  
 compress\_zip() (in module aquaduct.utils.helpers), 69  
 concatenate() (in module aquaduct.utils.helpers), 69  
 ConnectToPymol (class in aquaduct.visual.pymol\_connector), 72  
 containing\_residues() (AtomSelection method), 61  
 contains() (SmartRangeFunction method), 69  
 contains\_residues() (AtomSelection method), 61  
 convert() (LoadDumpWrapper method), 37  
 convexhull() (SinglePathPlotter method), 73  
 coords (GenericPaths attribute), 56  
 coords (Inlets attribute), 54  
 coords (MacroMolPath attribute), 57  
 coords (PassingPath attribute), 58  
 coords() (AtomSelection method), 61  
 coords() (FakeSingleResidueSelection method), 43  
 coords() (ResidueSelection method), 61  
 coords() (Selection method), 61  
 coords() (SingleResidueSelection method), 62  
 coords\_cont (MacroMolPath attribute), 57  
 coords\_filo (MacroMolPath attribute), 57  
 coords\_first\_in (MacroMolPath attribute), 57  
 coords\_first\_in (PassingPath attribute), 58  
 coords\_in (MacroMolPath attribute), 57  
 coords\_last\_out (MacroMolPath attribute), 57  
 coords\_last\_out (PassingPath attribute), 58  
 coords\_object (MacroMolPath attribute), 57  
 coords\_out (MacroMolPath attribute), 57  
 coords\_range() (in module aquaduct.traj.sandwich), 62  
 coords\_range\_core() (in module aquaduct.traj.sandwich), 62  
 coords\_smooth() (FakeSingleResidueSelection method), 43  
 coords\_smooth() (SingleResidueSelection method), 62  
 coords\_types\_prob\_widths() (CTypeSpathsCollectionWorker method), 41  
 CoordsRangeIndexCache (class in aquaduct.apps.data), 37  
 correct\_window() (MasterReader method), 59

cpu\_count (CpuThreadsCount attribute), 71  
CpuThreadsCount (class in aquaduct.utils.multip), 71  
create\_tmpfile() (in module aquaduct.utils.helpers), 66  
ct\_file (ConnectToPymol attribute), 72  
ct\_pymol (ConnectToPymol attribute), 72  
CTypeSpathsCollection (class in aquaduct.geom.master), 41  
CTypeSpathsCollectionWorker (class in aquaduct.geom.master), 41  
cut\_thyself() (WhereToCut method), 53

## D

decode\_color() (ConnectToPymol method), 73  
del\_cgo\_object() (ConnectToPymol method), 73  
derivative() (in module aquaduct.geom.traces), 48  
dict2arrays() (IdsOverIds static method), 38  
diff() (in module aquaduct.geom.traces), 48  
distance() (ColorMapDistMap method), 71  
DistanceWindowOverMaxStepSmooth (class in aquaduct.geom.smooth), 47  
DistanceWindowSmooth (class in aquaduct.geom.smooth), 46  
dump() (ValveDataAccess\_pickle method), 38  
dump\_config() (ValveConfig method), 39  
dump\_frames() (ReaderTraj method), 60  
dump\_frames() (TmpDumpWriterOfMDA method), 53  
dump\_frames\_to\_file() (ReaderTraj method), 60  
dump\_frames\_to\_file() (ReaderTrajViaMDA method), 61

## E

emit\_message\_to\_file\_in\_root\_logger() (in module aquaduct.utils.clui), 62  
empty\_coords (MacroMolPath attribute), 56  
ends (MacroMolPath attribute), 57  
engine (MasterReader attribute), 59  
engine\_name (MasterReader attribute), 59  
ETA() (SimpleProgressBar method), 64  
etypes (MacroMolPath attribute), 57  
etypes\_cont (MacroMolPath attribute), 57  
euclidean() (in module aquaduct.visual.helpers), 71  
extend\_inlets() (Inlets method), 54

## F

f\_like() (in module aquaduct.visual.helpers), 71  
FakeSingleResidueSelection (class in aquaduct.geom.master), 43  
fbm (class in aquaduct.utils.clui), 62  
find\_paths() (GenericPaths method), 56  
find\_paths\_types() (GenericPaths method), 56  
finish() (SimpleProgressBar method), 65  
first\_element() (SmartRange method), 70  
first\_element() (SmartRangeFunction method), 69  
fit() (BarberCluster method), 40  
fit() (PerformClustering method), 40  
float\_default (NumpyDefaultsStorageTypes attribute), 71  
FloatWindow (class in aquaduct.geom.smooth), 45

frames (GenericPaths attribute), 56  
FramesRangeCollection (class in aquaduct.traj.sandwich), 62  
full\_size() (CTypeSpathsCollection method), 42

## G

GeneralWindow (class in aquaduct.geom.smooth), 45  
generic (InletClusterExtendedType attribute), 54  
GenericPaths (class in aquaduct.traj.paths), 56  
GenericPathTypeCodes (class in aquaduct.traj.paths), 55  
get() (BasicPymolCGO method), 72  
get() (SmartRange method), 70  
get() (SmartRangeDecrement method), 70  
get() (SmartRangeEqual method), 69  
get() (SmartRangeFunction method), 69  
get() (SmartRangeIncrement method), 70  
get\_acceleration\_cont() (MacroMolPath method), 58  
get\_ax3d() (in module aquaduct.visual.quickplot), 73  
get\_chull() (Inlets method), 55  
get\_cluster\_options() (ValveConfig method), 39  
get\_cmap() (in module aquaduct.visual.helpers), 71  
get\_common\_traj\_data() (ValveConfig method), 39  
get\_coords() (MacroMolPath method), 58  
get\_coords() (PassingPath method), 58  
get\_coords\_cont() (MacroMolPath method), 58  
get\_cric\_reader() (in module aquaduct.apps.data), 37  
get\_current\_nr() (WhereToCut method), 53  
get\_default\_config() (ValveConfig method), 39  
get\_distance\_both\_cont() (MacroMolPath method), 58  
get\_distance\_cont() (MacroMolPath method), 58  
get\_distance\_rev\_cont() (MacroMolPath method), 58  
get\_flat\_tree() (Inlets method), 55  
get\_general\_comment() (ValveConfig method), 39  
get\_global\_options() (ValveConfig method), 39  
get\_inlets() (MacroMolPath method), 57  
get\_inlets() (PassingPath method), 58  
get\_inlets\_references() (Inlets method), 54  
get\_leaf() (SimpleTree method), 65  
get\_master\_path() (CTypeSpathsCollection method), 43  
get\_object\_from\_name() (in module aquaduct.apps.data), 38  
get\_object\_name() (in module aquaduct.apps.data), 38  
get\_paths\_in() (GenericPaths method), 56  
get\_paths\_out() (GenericPaths method), 56  
get\_ranges() (FramesRangeCollection method), 62  
get\_reader() (Selection method), 61  
get\_reader() (SingleResidueSelection method), 62  
get\_reader\_by\_id() (MasterReader method), 59  
get\_real() (Window method), 59  
get\_recluster\_options() (ValveConfig method), 39  
get\_required\_params() (in module aquaduct.geom.cluster), 39  
get\_single\_path\_types() (GenericPaths method), 56  
get\_single\_reader() (MasterReader method), 59  
get\_smooth\_options() (ValveConfig method), 39  
get\_stage\_options() (ValveConfig method), 39



- get\_str\_timestamp() (in module aquaduct.utils.clui), 65
  - get\_variable() (ValveDataAccess\_pickle method), 38
  - get\_vda\_reader() (in module aquaduct.apps.data), 37
  - get\_velocity\_cont() (MacroMolPath method), 58
  - global\_name() (ValveConfig static method), 38
  - GlobalConfigStore (class in aquaduct.apps.data), 37
  - glue() (in module aquaduct.traj.paths), 55
  - greetings() (in module aquaduct), 74
  - gregorian\_year\_in\_days (in module aquaduct.utils.clui), 63
  - grey (ColorMapDistMap attribute), 71
  - gsep() (in module aquaduct.utils.clui), 63
  - gtypes (MacroMolPath attribute), 57
  - gtypes (PassingPath attribute), 58
  - gtypes\_cont (MacroMolPath attribute), 57
- ## H
- has\_in (MacroMolPath attribute), 57
  - has\_in (PassingPath attribute), 58
  - has\_object (MacroMolPath attribute), 58
  - has\_out (MacroMolPath attribute), 58
  - has\_out (PassingPath attribute), 58
  - heartbeat() (SimpleProgressBar method), 65
  - here() (LinearizeOneWay method), 50
  - here() (LinearizeRecursive method), 50
- ## I
- ids() (Selection method), 61
  - IdsOverIds (class in aquaduct.apps.data), 38
  - incoming (ProtoInletTypeCodes attribute), 53
  - init\_ax() (MPLTracePlotter method), 74
  - init\_pymol() (ConnectToPymol method), 73
  - init\_script() (ConnectToPymol method), 73
  - Inlet (class in aquaduct.traj.inlets), 54
  - inlet2sphere() (WhereToCut method), 53
  - InletClusterExtendedType (class in aquaduct.traj.inlets), 54
  - InletClusterGenericType (class in aquaduct.traj.inlets), 53
  - Inlets (class in aquaduct.traj.inlets), 54
  - InletTypeCodes (class in aquaduct.traj.inlets), 53
  - input (InletClusterGenericType attribute), 54
  - int2range() (in module aquaduct.utils.helpers), 67
  - int\_default (NumpyDefaultsStorageTypes attribute), 71
  - internal (ProtoInletTypeCodes attribute), 53
  - internal\_incoming (InletTypeCodes attribute), 53
  - internal\_outgoing (InletTypeCodes attribute), 53
  - IntWindow (class in aquaduct.geom.smooth), 45
  - is\_frame\_in() (MacroMolPath method), 57
  - is\_frame\_object() (MacroMolPath method), 57
  - is\_frame\_out() (MacroMolPath method), 57
  - is\_frame\_walk() (MacroMolPath method), 57
  - is\_frame\_walk() (PassingPath method), 58
  - is\_iterable() (in module aquaduct.utils.helpers), 67
  - is\_leaf() (SimpleTree method), 65
  - is\_linear() (TriangleLinearize method), 50
  - is\_linear() (VectorLinearize method), 51
  - is\_linear\_core() (VectorLinearize method), 51
  - is\_number() (in module aquaduct.utils.helpers), 66
  - is\_overlapping\_with\_cloud() (WhereToCut method), 53
  - is\_passing() (MacroMolPath method), 57
  - is\_passing() (PassingPath method), 58
  - is\_passing() (SinglePath method), 58
  - is\_point\_within() (Sphere method), 52
  - is\_point\_within\_convexhull() (in module aquaduct.geom.convexhull), 40
  - is\_single() (MacroMolPath method), 57
  - is\_single() (PassingPath method), 58
  - is\_single() (SinglePath method), 58
  - is\_sphere\_cloud() (Sphere method), 52
  - is\_sphere\_within() (Sphere method), 52
  - isin() (SmartRange method), 70
  - isin() (SmartRangeDecrement method), 70
  - isin() (SmartRangeEqual method), 70
  - isin() (SmartRangeFunction method), 69
  - isin() (SmartRangeIncrement method), 70
  - iterate() (MasterReader method), 59
  - iterate\_or\_die() (in module aquaduct.utils.helpers), 69
  - iterate\_over\_frames() (ReaderTraj method), 60
  - itype (InletTypeCodes attribute), 53
  - ix() (Selection method), 61
- ## L
- last\_element() (SmartRange method), 70
  - last\_element() (SmartRangeDecrement method), 70
  - last\_element() (SmartRangeEqual method), 70
  - last\_element() (SmartRangeFunction method), 69
  - last\_element() (SmartRangeIncrement method), 70
  - last\_times() (SmartRange method), 70
  - layer() (Selection method), 61
  - leafs\_names (SimpleTree attribute), 65
  - left() (in module aquaduct.traj.paths), 55
  - len() (Selection method), 61
  - len() (Window method), 59
  - length\_step\_std() (in module aquaduct.geom.traces), 48
  - lens() (CTypeSpathsCollection method), 42
  - lens\_norm() (CTypeSpathsCollection method), 42
  - lens\_real() (CTypeSpathsCollection method), 42
  - lim2clusters() (Inlets method), 55
  - lim2rnames() (Inlets method), 55
  - lim2spaths() (Inlets method), 55
  - lim2types() (Inlets method), 55
  - lim\_to() (Inlets method), 55
  - limspaths2() (Inlets method), 55
  - lind() (in module aquaduct.utils.helpers), 66
  - LinearizeHobbit (class in aquaduct.geom.traces), 50
  - LinearizeHobbitTriangle (class in aquaduct.geom.traces), 51
  - LinearizeHobbitVector (class in aquaduct.geom.traces), 51
  - LinearizeOneWay (class in aquaduct.geom.traces), 50
  - LinearizeOneWayTriangle (class in aquaduct.geom.traces), 52
  - LinearizeOneWayVector (class in aquaduct.geom.traces), 51
  - LinearizeRecursive (class in aquaduct.geom.traces), 50

LinearizeRecursiveTriangle (class in aquaduct.geom.traces), 51  
LinearizeRecursiveVector (class in aquaduct.geom.traces), 51  
list\_blocks\_to\_slices() (in module aquaduct.utils.helpers), 68  
listify() (in module aquaduct.utils.helpers), 68  
load() (ValveDataAccess\_pickle method), 38  
load\_config() (ValveConfig method), 39  
load\_cric() (in module aquaduct.apps.data), 37  
load\_pdb() (ConnectToPymol method), 73  
LoadDumpWrapper (class in aquaduct.apps.data), 37

## M

MacroMolPath (class in aquaduct.traj.paths), 56  
make\_color\_triple() (BasicPymolCGO static method), 72  
make\_default\_array() (in module aquaduct.utils.maths), 71  
make\_iterable() (in module aquaduct.utils.helpers), 69  
make\_val() (InletClusterGenericType method), 54  
MasterPath (class in aquaduct.traj.paths), 59  
MasterReader (class in aquaduct.traj.sandwich), 59  
max() (SmartRange method), 70  
max\_frame (GenericPaths attribute), 56  
max\_window\_at\_pos() (GeneralWindow static method), 45  
MaxStepSmooth (class in aquaduct.geom.smooth), 46  
MeanShiftBandwidth() (in module aquaduct.geom.cluster), 40  
message() (in module aquaduct.utils.clui), 62  
message\_special() (in module aquaduct.utils.clui), 62  
midpoints() (in module aquaduct.geom.traces), 48  
mimic\_old\_var\_name (ValveDataAccess\_pickle attribute), 38  
min() (SmartRange method), 70  
min\_frame (GenericPaths attribute), 56  
MPLTracePlotter (class in aquaduct.visual.quickplot), 74

## N

names() (ResidueSelection method), 61  
new() (BasicPymolCGO method), 72  
next() (SimpleProgressBar method), 65  
noaction() (in module aquaduct.utils.helpers), 68  
Normalize (class in aquaduct.geom.pca), 43  
number\_of\_frames() (MasterReader method), 60  
number\_of\_layers() (MasterReader method), 60  
numbers() (Selection method), 61  
NumpyDefaultsStorageTypes (class in aquaduct.utils.maths), 71

## O

object\_len (MacroMolPath attribute), 57  
object\_len (PassingPath attribute), 58  
object\_name (GenericPathTypeCodes attribute), 56  
open() (SimpleTarWriteHelper method), 72  
open() (ValveDataAccess\_nc method), 38

open() (ValveDataAccess\_pickle method), 38  
open() (ValveDataAccessRoots method), 38  
open\_reader\_traj (MasterReader attribute), 59  
open\_trajectory() (ReaderTraj method), 60  
open\_trajectory() (ReaderTrajViaMDA method), 60  
orient\_on() (ConnectToPymol method), 73  
out\_name (GenericPathTypeCodes attribute), 56  
outgoing (ProtoInletTypeCodes attribute), 53  
output (InletClusterGenericType attribute), 54  
overlaps() (SmartRangeFunction method), 69  
overlaps\_mutual() (SmartRangeFunction method), 69

## P

P (PCA attribute), 44  
parse\_selection() (ReaderTraj method), 60  
parse\_selection() (ReaderTrajViaMDA method), 61  
part2type\_dict (in module aquaduct.geom.master), 41  
parts (CTypeSpathsCollection attribute), 42  
parts (in module aquaduct.geom.master), 41  
PassingPath (class in aquaduct.traj.paths), 58  
path (PassingPath attribute), 58  
path\_in (MacroMolPath attribute), 57  
path\_in\_code (PathTypesCodes attribute), 55  
path\_object (MacroMolPath attribute), 57  
path\_object\_code (PathTypesCodes attribute), 55  
path\_out (MacroMolPath attribute), 57  
path\_out\_code (PathTypesCodes attribute), 55  
path\_trace() (SimpleTracePlotter method), 73  
path\_walk\_code (PathTypesCodes attribute), 55  
paths (MacroMolPath attribute), 57  
paths (PassingPath attribute), 58  
paths\_cont (MacroMolPath attribute), 57  
paths\_first\_in (MacroMolPath attribute), 57  
paths\_first\_in (PassingPath attribute), 58  
paths\_inlets() (SinglePathPlotter method), 73  
paths\_last\_out (MacroMolPath attribute), 57  
paths\_last\_out (PassingPath attribute), 58  
paths\_trace() (SinglePathPlotter method), 73  
PathTypesCodes (class in aquaduct.traj.paths), 55  
pbar (in module aquaduct.utils.clui), 65  
PCA (class in aquaduct.geom.pca), 44  
percent() (SimpleProgressBar method), 64  
perform\_clustering() (Inlets method), 55  
perform\_reclustering() (Inlets method), 55  
PerformClustering (class in aquaduct.geom.cluster), 40  
plot\_colorful\_lines() (in module aquaduct.visual.quickplot), 73  
plot\_line() (MPLTracePlotter method), 74  
plot\_line() (SimpleTracePlotter method), 73  
plot\_spath\_spectrum() (in module aquaduct.visual.quickplot), 73  
preprocess() (PCA method), 44  
preprocess\_undo() (PCA method), 44  
print\_simple\_tree() (in module aquaduct.utils.clui), 65  
protein\_trace() (SimpleProteinPlotter method), 73  
ProtoInletTypeCodes (class in aquaduct.traj.inlets), 53



## R

[range\(\)](#) (Window method), 59  
[range2int\(\)](#) (in module aquaduct.utils.helpers), 67  
[raw](#) (SmartRange attribute), 70  
[read\(\)](#) (LoadDumpWrapper method), 37  
[reader](#) (ReaderAccess attribute), 60  
[ReaderAccess](#) (class in aquaduct.traj.sandwich), 60  
[ReaderTraj](#) (class in aquaduct.traj.sandwich), 60  
[ReaderTrajViaMDA](#) (class in aquaduct.traj.sandwich), 60  
[readline\(\)](#) (LoadDumpWrapper method), 37  
[real\\_number\\_of\\_frames\(\)](#) (MasterReader method), 60  
[real\\_number\\_of\\_frames\(\)](#) (ReaderTraj method), 60  
[real\\_number\\_of\\_frames\(\)](#) (ReaderTrajViaMDA method), 61  
[recluster\\_cluster\(\)](#) (Inlets method), 55  
[recluster\\_name\(\)](#) (ValveConfig static method), 38  
[recluster\\_outliers\(\)](#) (Inlets method), 55  
[recursive\\_clusterization\\_name\(\)](#) (ValveConfig static method), 38  
[recursive\\_threshold\\_name\(\)](#) (ValveConfig static method), 38  
[refs](#) (Inlets attribute), 54  
[refs\\_names](#) (Inlets attribute), 54  
[renumber\\_clusters\(\)](#) (Inlets method), 55  
[residues\(\)](#) (AtomSelection method), 61  
[residues\\_names\(\)](#) (ReaderTraj method), 60  
[residues\\_names\(\)](#) (ReaderTrajViaMDA method), 61  
[residues\\_positions\(\)](#) (ReaderTraj method), 60  
[residues\\_positions\(\)](#) (ReaderTrajViaMDA method), 61  
[ResidueSelection](#) (class in aquaduct.traj.sandwich), 61  
[resize\\_leaf\\_0\(\)](#) (Inlets method), 54  
[rev\(\)](#) (SmartRange method), 70  
[rev\(\)](#) (SmartRangeDecrement method), 70  
[rev\(\)](#) (SmartRangeEqual method), 70  
[rev\(\)](#) (SmartRangeFunction method), 69  
[rev\(\)](#) (SmartRangeIncrement method), 70  
[right\(\)](#) (in module aquaduct.traj.paths), 55  
[robust\\_and\(\)](#) (in module aquaduct.utils.helpers), 66  
[robust\\_or\(\)](#) (in module aquaduct.utils.helpers), 66  
[roman\\_emulation](#) (class in aquaduct.utils.clui), 62  
[roots](#) (ValveDataAccessRoots attribute), 38  
[rotate](#) (SimpleProgressBar attribute), 64

## S

[sandwich\(\)](#) (MasterReader method), 59  
[sandwich\\_mode](#) (MasterReader attribute), 59  
[save\\_config\(\)](#) (ValveConfig method), 39  
[save\\_config\\_stream\(\)](#) (ValveConfig method), 39  
[save\\_cric\(\)](#) (in module aquaduct.apps.data), 37  
[save\\_file2tar\(\)](#) (SimpleTarWriteHelper method), 72  
[save\\_object2tar\(\)](#) (SimpleTarWriteHelper method), 72  
[SavgolSmooth](#) (class in aquaduct.geom.smooth), 47  
[scatter\(\)](#) (MPLTracePlotter method), 74  
[scatter\(\)](#) (SinglePathPlotter method), 73  
[scope\\_name](#) (GenericPathTypeCodes attribute), 56  
[Selection](#) (class in aquaduct.traj.sandwich), 61  
[set\\_frame\(\)](#) (ReaderTraj method), 60

[set\\_real\\_frame\(\)](#) (ReaderTraj method), 60  
[set\\_real\\_frame\(\)](#) (ReaderTrajViaMDA method), 61  
[set\\_savgol\\_function\(\)](#) (SavgolSmooth method), 47  
[set\\_variable\(\)](#) (ValveDataAccess\_pickle method), 38  
[show\(\)](#) (SimpleProgressBar method), 65  
[showit\(\)](#) (in module aquaduct.visual.quickplot), 73  
[simple\\_types\\_distribution\(\)](#) (CTypeSpathsCollection static method), 42  
[SimplePathPlotter](#) (class in aquaduct.visual.quickplot), 74  
[SimpleProgressBar](#) (class in aquaduct.utils.clui), 64  
[SimpleProteinPlotter](#) (class in aquaduct.visual.quickplot), 73  
[SimpleTarWriteHelper](#) (class in aquaduct.visual.pymol\_connector), 72  
[SimpleTracePlotter](#) (class in aquaduct.visual.quickplot), 73  
[SimpleTree](#) (class in aquaduct.utils.clui), 65  
[single\\_path\\_traces\(\)](#) (SimplePathPlotter method), 74  
[single\\_residues\(\)](#) (ResidueSelection method), 61  
[single\\_trace\(\)](#) (SimpleTracePlotter method), 73  
[SinglePath](#) (class in aquaduct.traj.paths), 58  
[SinglePathID](#) (class in aquaduct.traj.paths), 56  
[SinglePathPlotter](#) (class in aquaduct.visual.pymol\_connector), 73  
[SingleResidueSelection](#) (class in aquaduct.traj.sandwich), 62  
[size](#) (Inlets attribute), 54  
[size](#) (MacroMolPath attribute), 57  
[sizes](#) (MacroMolPath attribute), 57  
[sizes](#) (PassingPath attribute), 58  
[small\\_clusters\\_to\\_outliers\(\)](#) (Inlets method), 55  
[smart\\_time\\_string\(\)](#) (in module aquaduct.utils.clui), 63  
[SmartRange](#) (class in aquaduct.utils.helpers), 70  
[SmartRangeDecrement](#) (class in aquaduct.utils.helpers), 70  
[SmartRangeEqual](#) (class in aquaduct.utils.helpers), 69  
[SmartRangeFunction](#) (class in aquaduct.utils.helpers), 69  
[SmartRangeIncrement](#) (class in aquaduct.utils.helpers), 70  
[Smooth](#) (class in aquaduct.geom.smooth), 44  
[smooth\(\)](#) (ActiveWindowOverMaxStepSmooth method), 47  
[smooth\(\)](#) (ActiveWindowSmooth method), 46  
[smooth\(\)](#) (DistanceWindowOverMaxStepSmooth method), 48  
[smooth\(\)](#) (DistanceWindowSmooth method), 46  
[smooth\(\)](#) (MaxStepSmooth method), 47  
[smooth\(\)](#) (SavgolSmooth method), 47  
[smooth\(\)](#) (Smooth method), 44  
[smooth\(\)](#) (WindowOverMaxStepSmooth method), 47  
[smooth\(\)](#) (WindowSmooth method), 46  
[smooth\\_coords\\_ranges\(\)](#) (in module aquaduct.traj.sandwich), 62  
[smooth\\_name\(\)](#) (ValveConfig static method), 39  
[sort\\_clusters\(\)](#) (Inlets method), 55  
[sortify\(\)](#) (in module aquaduct.utils.helpers), 67

`spath2ctype()` (Inlets method), 55  
`spath2spheres()` (WhereToCut method), 53  
`spath_spectrum()` (in module `aquaduct.visual.quickplot`), 73  
`spaths2ctypes()` (Inlets method), 55  
`spaths_spectra()` (in module `aquaduct.visual.quickplot`), 73  
`Sphere` (class in `aquaduct.traj.barber`), 52  
`split_list()` (in module `aquaduct.utils.helpers`), 69  
`stage_I_run()` (in module `aquaduct.apps.valvecore`), 39  
`stage_II_run()` (in module `aquaduct.apps.valvecore`), 39  
`stage_III_run()` (in module `aquaduct.apps.valvecore`), 39  
`stage_IV_run()` (in module `aquaduct.apps.valvecore`), 39  
`stage_names()` (ValveConfig method), 39  
`stage_V_run()` (in module `aquaduct.apps.valvecore`), 39  
`stage_VI_run()` (in module `aquaduct.apps.valvecore`), 39  
`Standartize` (class in `aquaduct.geom.pca`), 44  
`strech_zip()` (in module `aquaduct.utils.helpers`), 69  
`surface` (ProtoInletTypeCodes attribute), 53  
`surface_incoming` (InletTypeCodes attribute), 53  
`surface_outgoing` (InletTypeCodes attribute), 53

## T

`thead()` (in module `aquaduct.utils.clui`), 64  
`threads_count` (CpuThreadsCount attribute), 71  
`tictoc` (class in `aquaduct.utils.clui`), 63  
`TmpDumpWriterOfMDA` (class in `aquaduct.traj.dumps`), 53  
`topology` (MasterReader attribute), 59  
`topology_resid` (SingleResidueSelection attribute), 62  
`topology_resids()` (ReaderTraj method), 60  
`topology_resids()` (ReaderTrajViaMDA method), 61  
`toRoman()` (roman\_emulation method), 62  
`tracepoints()` (in module `aquaduct.geom.traces`), 48  
`trajectory` (MasterReader attribute), 59  
`triangle_angles()` (in module `aquaduct.geom.traces`), 48  
`triangle_angles_last()` (in module `aquaduct.geom.traces`), 48  
`triangle_height()` (in module `aquaduct.geom.traces`), 49  
`TriangleLinearize` (class in `aquaduct.geom.traces`), 50  
`tsep()` (in module `aquaduct.utils.clui`), 63  
`ttime()` (SimpleProgressBar method), 65  
`tupleify()` (in module `aquaduct.utils.helpers`), 68  
`type` (SmartRangeDecrement attribute), 70  
`type` (SmartRangeEqual attribute), 69  
`type` (SmartRangeIncrement attribute), 70  
`types` (GenericPaths attribute), 56  
`types` (Inlets attribute), 54  
`types` (MacroMolPath attribute), 57  
`types` (PassingPath attribute), 58  
`types_cont` (MacroMolPath attribute), 57  
`types_distribution()` (CTypeSpathsCollection method), 43  
`types_in` (MacroMolPath attribute), 57  
`types_object` (MacroMolPath attribute), 57

`types_out` (MacroMolPath attribute), 57  
`types_prob_to_types()` (CTypeSpathsCollection method), 43

## U

`underline()` (in module `aquaduct.utils.clui`), 64  
`undo()` (Center method), 43  
`undo()` (Normalize method), 44  
`undo()` (PCA method), 44  
`undo()` (Standartize method), 44  
`union()` (in module `aquaduct.traj.paths`), 55  
`union_full()` (in module `aquaduct.traj.paths`), 55  
`union_smartr()` (in module `aquaduct.traj.paths`), 55  
`uniqify()` (in module `aquaduct.utils.helpers`), 68  
`uniqify()` (Selection method), 61  
`unknown_names` (ValveDataAccess\_pickle attribute), 38  
`update()` (CTypeSpathsCollection method), 42  
`update()` (SimpleProgressBar method), 65

## V

`valve_begin()` (in module `aquaduct.apps.valvecore`), 39  
`valve_end()` (in module `aquaduct.apps.valvecore`), 39  
`valve_exec_stage()` (in module `aquaduct.apps.valvecore`), 39  
`valve_load_config()` (in module `aquaduct.apps.valvecore`), 39  
`ValveConfig` (class in `aquaduct.apps.valvecore`), 38  
`ValveDataAccess` (in module `aquaduct.apps.data`), 38  
`ValveDataAccess_nc` (class in `aquaduct.apps.data`), 38  
`ValveDataAccess_pickle` (class in `aquaduct.apps.data`), 37  
`ValveDataAccessRoots` (class in `aquaduct.apps.data`), 38  
`vdw()` (AtomSelection method), 61  
`VdW_radii` (in module `aquaduct.traj.sandwich`), 60  
`vector_norm()` (in module `aquaduct.geom.traces`), 48  
`VectorLinearize` (class in `aquaduct.geom.traces`), 50  
`vectors_angle()` (in module `aquaduct.geom.traces`), 49  
`vectors_angle_alt()` (in module `aquaduct.geom.traces`), 49  
`vectors_angle_alt_anorm()` (in module `aquaduct.geom.traces`), 49  
`vectors_angle_anorm()` (in module `aquaduct.geom.traces`), 49  
`version()` (in module `aquaduct`), 74  
`version_nice()` (in module `aquaduct`), 74

## W

`what2what()` (in module `aquaduct.utils.helpers`), 69  
`WhereToCut` (class in `aquaduct.traj.barber`), 52  
`Window` (class in `aquaduct.traj.sandwich`), 59  
`window` (MasterReader attribute), 59  
`WindowOverMaxStepSmooth` (class in `aquaduct.geom.smooth`), 47  
`WindowSmooth` (class in `aquaduct.geom.smooth`), 46

## X

`xor()` (in module `aquaduct.traj.paths`), [55](#)  
`xor_full()` (in module `aquaduct.traj.paths`), [55](#)  
`xor_smartr()` (in module `aquaduct.traj.paths`), [55](#)  
`xzip_xzip()` (in module `aquaduct.utils.helpers`), [69](#)

## Y

`yield_generic_paths()` (in module `aquaduct.traj.paths`),  
[56](#)  
`yield_single_paths()` (in module `aquaduct.traj.paths`),  
[56](#)  
`yield_spath_len_and_smooth_diff_in_types_slices()`  
(in module `aquaduct.visual.quickplot`), [73](#)

## Z

`zip_zip()` (in module `aquaduct.utils.helpers`), [69](#)